

PHILIPS

InnoSpire Go

InnoSpire Go

Mesh nebulizer system

Index

| | | |
|-------------------|----|-----|
| English | EN | 4 |
| ภาษาไทย | TH | 15 |
| العربية | AR | 26 |
| תירבע | IL | 37 |
| Bahasa Melayu | MS | 47 |
| Dansk | DA | 57 |
| Deutsch | DE | 67 |
| Nederlands | NL | 78 |
| Français | FR | 88 |
| Bahasa Indonesian | ID | 99 |
| Italiano | IT | 109 |
| Norsk | NO | 120 |
| Português | PT | 130 |
| Español | ES | 141 |
| Suomi | FI | 152 |
| Svenska | SV | 163 |
| Русский | RU | 173 |
| Ελληνικά | EL | 184 |
| EMC data tables | EN | 196 |

- | | |
|--|---------------------------------------|
| 1. Handset | 7. Medication chamber lid |
| 2. Mouthpiece assembly | 8. Medication chamber |
| 3. AC power adapter | 9. Mouthpiece assembly release button |
| 4. LiteTouch medium mask (age 1 – 5 years) | 10. Power socket (back view) |
| 5. Mask adapter | 11. Carry case (not shown) |
| 6. On/Off button and LED indicator | |

Instructions for use

Please read these instructions carefully before first use. If you do not understand any portion of these instructions, contact your healthcare provider or call Philips Customer Service at 1-724-387-4000.

General information

Intended use: The InnoSpire Go is a general purpose mesh nebulizer that is intended to be used to nebulize commonly prescribed liquid inhaled medications for respiratory disease. It is intended for single patient use, to deliver multiple doses.

It is to be used by patients from infant to adult.

It is suitable for use in the home environment, but may also be used in a hospital/clinic setting.

 **Keep these Instructions for future reference**

Retain carton and packing materials for storing the unit or for product returns.

Cautions

It is recommended to have a backup device (e.g. MDI) for respiratory delivery in case a situation arises when your nebulizer cannot be used.


- It is not to be used by patients that are unconscious or not breathing spontaneously.
- The patient or the patient's care giver is the intended operator of the device.
- Only use this nebulizer with medications prescribed by your physician.
- Do not place or store the product where it can fall into water.
- Do not submerge the handset in water or other liquid cleaning agents.
- If the handset is accidentally dropped into a liquid, immediately disconnect the AC adapter from the wall outlet, prior to removing the handset from the liquid.
- Any liquid spilled on the handset should be allowed to dry before operating.
- Do not position the product where it is difficult to disconnect the plug.
- Adult supervision is necessary when this product is used by, on, or near children or physically challenged individuals or individuals with learning difficulties.
- Use this product only for its intended use as described in this manual. Failure to do so may compromise performance. Only use accessories provided with the nebulizer and recommended by the manufacturer. Store in a clean place out of the reach of children.

- Do not disassemble or modify the nebulizer in anyway. There are no serviceable parts. The battery is non-replaceable.
- Never charge or operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water.
- Do not connect to other equipment not described in these instructions.
- Do not charge your device on an aircraft.
- Keep the power cord away from heated surfaces.
- The nebulizer must be operated using the specified power sources.
- Always check the inside of the mouthpiece assembly or mask for any debris before use.
- Do not use while operating a vehicle.
- When operated at an ambient temperature 40°C the handset can reach temperatures of up to 43°C. The device should not be used for more than 10 minutes during this scenario.
- Do not use in an anesthetic or ventilator breathing system.
- InnoSpire Go and all its parts (including battery) must be disposed of properly and according to the local regulations in force (for example the WEEE directive).
- Small parts can be inhaled or swallowed. In addition, the cable, due to its length, may result in provoking strangulation or asphyxiation. Do not leave the device alone with a small child or physically challenged individual or individual with learning difficulties.
- The nebulizer should not be used around flammable substances e.g. oxygen, nitrous oxide, or in the presence of a flammable anesthetic mixture.
- Never poke the mesh, or clean it with any sharp objects. This could damage the mesh and prevent your device from operating properly.
- Do not autoclave the device.
- Report unexpected operation or events to Philips.
- Precautions are to be taken in the event of changes in the performance of the device, please refer to the troubleshooting section.

How to use your InnoSpire Go

- 1 After unpacking the nebulizer check you have all the items listed and that there is no visible damage or defects. Contact your product distributor or Philips customer service if anything is missing or damaged. Prior to first use clean the mouthpiece assembly as per the cleaning instructions and fully charge the battery. Ensure the device is disconnected from the power during assembly and disassembly.
- 2 Check the nebulizer and accessories are clean, dry and not damaged prior to use.
- 3 Attach the mouthpiece assembly to the handset. **Do not use separate mask adapter as a mouthpiece.**
- 4 Lift the green medication chamber lid.

- 5 Empty the contents of the medication vial into the medication chamber.
 - The maximum fill volume is 8ml and this is shown by the word MAX on the protrusion below the hinge. **Do not fill medication above this level.**
- 6 Close the lid to the medication chamber.
- 7 If using a mask, attach the mask adapter to the mouthpiece assembly and attach the mask to mask adapter.
- 8 If you are using the power adapter to power the nebulizer, plug the cable into the socket on the handset and plug the adapter into the wall outlet.
- 9 Press the on/off button on the handset to switch the nebulizer on and begin nebulization.
- 10 Check the battery level.
 - If the LED is **SOLID GREEN** the battery is charged.
 - If the LED is **SOLID AMBER**, there is enough charge for at least one more treatment. Please charge your nebulizer after your treatment.
 - If the LED **FLASHES AMBER** and then switches off there is not enough charge to take a treatment.
 - If the LED does not illuminate, please refer to the troubleshooting section of this manual.
- 11 Ensure aerosol is coming out of the mouthpiece or mask.
- 12 a) **If using the mouthpiece**, hold the handset in your hand and place the mouthpiece between your teeth with your lips sealed around it. Breathe normally through your mouth.
 b) **If using a mask**, hold the handset in your hand and gently press the mask against your face and breathe normally through your mouth.
 During use, some aerosol will be emitted from the back of the mouthpiece assembly.


Do not tilt the device in any direction more than 45 degrees during the course of a treatment as this will prevent the nebulizer from completely nebulizing all the medication in the chamber.

- 13 If you need to take a rest, press the on/off button to stop your treatment. To continue your treatment press the on/off button again.
- 14 Your treatment is finished when the nebulizer beeps and the LED flashes. The device will turn off automatically.
- 15 Check the medication chamber for residual medication. If there are more than a few drops remaining, press the on/off button again to continue your treatment.
- 16 Clean the nebulizer following the cleaning instructions.

Charging the battery

- 1 Connect the power adapter to the socket on the back of the handset.
- 2 Plug the other end of the power adapter into the wall outlet.
- 3 The LED will **PULSE GREEN**.
- 4 Charge the battery until the LED turns **SOLID GREEN** indicating it is fully charged.
- 5 Unplug the power adapter from both the handset and the wall outlet. It is recommended to unplug the power adapter once the battery is fully charged to preserve battery life.

Cleaning and maintenance

⚠ Cautions

- Do not autoclave the mouthpiece assembly or handset.
- Do not put any part of the device in a microwave or conventional oven.
- Do not immerse the handset in liquid or steam clean.
- Never poke or clean the mesh with any sharp objects as this will damage the mesh.
- Do not clean the device when in use.
- Disconnect device from power supply prior to cleaning.

Cleaning

After each use:

- Pour away residual medication from the medication chamber.
- Press the mouthpiece assembly release button to separate the mouthpiece assembly from the handset.
- Rinse the mouthpiece assembly thoroughly under running tap water.
- Shake off excess water and allow to air dry fully before storing.


Daily cleaning

- Wash the mouthpiece assembly by hand in a bowl of warm soapy water (liquid dishwashing soap) for 2 minutes.
- Rinse the mouthpiece assembly thoroughly under running tap water.
- Shake off excess water and allow to air dry fully before storing.

Weekly disinfection

CAUTION: Risk of scalding. Use care around boiling water and in handling hot parts.

- Prior to disinfection, ensure all parts are visibly clean and free from dirt.
- Boil the mouthpiece assembly in water for 10 minutes. Ensure the medication chamber lid is open and there is enough water in the pan to prevent the mouthpiece assembly from touching the bottom or boiling dry. Shake off excess water and allow to air dry fully before storing.

Or

- Disinfect the mouthpiece assembly in a Philips Avent baby bottle electric steam sterilizer, or equivalent, for one cycle (follow the instructions of the steam sterilizer). Allow to air dry fully before storing.

Or

- Wash the mouthpiece assembly by hand in a bowl of warm soapy water (liquid dishwashing soap) for 2 minutes. Rinse the mouthpiece assembly thoroughly under running tap water.
- Immerse the mouthpiece assembly in a disinfectant of the Gluteraldehyde group (testing performed with Korsolex Extra, 4% for 15 minutes), rinse thoroughly according to manufacturers instructions.
- Wipe your handset clean with a clean damp cloth weekly or as needed. Do not use other cleaning methods or solutions.

Mask cleaning

- If using a mask and mask adapter, once a week, agitate in warm soapy water for 2 minutes or soak for 10 minutes, rinse and air dry.

Maintenance

Mouthpiece assembly can be cleaned (boiled/disinfected) up to 52 times over 12 months.

To keep your InnoSpire Go working at the optimum level you should replace the mouthpiece assembly and mask, if used, every 12 months as this is a consumable part.

This is based on a typical usage profile of 3 treatments per day.

LED Indication

| LED Indication | What it means |
|-----------------------------------|--|
| SOLID GREEN | The battery is charged |
| SOLID AMBER | There is enough charge for at least one more treatment |
| FLASHES AMBER FIVE TIMES | There is not enough charge to take a treatment |
| FLASHES GREEN OR AMBER FIVE TIMES | The treatment is finished |
| PULSES GREEN | The battery is being charged |

Troubleshooting

| Problems | Solutions |
|---|--|
| The on button is pressed but nothing happens (no light, no nebulization). | Battery has insufficient charge, follow battery charging instructions. If the problem persists, contact your customer services representative. |
| LED illuminates on button press but goes off when button released. | Check to ensure that the mouthpiece assembly is attached properly. Make sure that the contacts on the handset are dry and free of debris. Make sure that there is medication loaded in the chamber. Recharge the battery. |
| When charging the battery, the LED does not pulse or turn on. | Ensure that the power cable is connected to the handset and plugged into the wall outlet or disconnect the power supply and reconnect. |
| Treatments take longer than usual with the same mouthpiece assembly. | Boil the mouthpiece assembly as per the weekly disinfection instructions. |
| Device indicates end of treatment but medication remains in medication chamber (more than a few drops). | Ensure nebulizer is held upright during nebulization. If the problem persists, contact your customer services representative. |
| Device does not indicate end of treatment even though all medication is nebulized. | Contact your customer services representative. |

Should the device still not operate properly after checking the unit as indicated above it may have a critical error, contact Philips Customer Service at 1-800-345-6443 or 1-724-387-4000.

Technical specifications

| | |
|--|--|
| Mains power supply | Input = 100 – 240 V~, 50/60 Hz Output = 5V \overleftrightarrow{DC} , 1.0A |
| Internal rechargeable battery power supply (lithium polymer) | 3.7 Volts nominal, 1200 mAh |
| Weight | 0.29 lbs/111g |
| Size | 7.0cm x 4.5cm x 13.5cm |

Class II device internally powered device (double safety insulation)
Type BF device (device with specific protection against electrical hazards)
Ingress Protection rating IP22. (Protected against solid foreign objects of 12.5mm diameter and greater; protected against vertically falling water drops when device titled upto 15°) The InnoSpire Go handset and charger form a Medical Electrical System, in which the charger is not Medical Electrical Equipment. The IP22 rating applies to the InnoSpire Go handset, the charger rating is IPX0.
All components shown in the illustration are applied parts.
Materials: Handset - Polyamide (PA) and Thermoplastic elastomer (TPE),
Mouthpiece Assembly - Polyamide (PA) and Polypropylene (PP),
Mask - polyester and silicone

Certification

| | |
|-------------------------------|---------------------------|
| Reference to standards | |
| Electric safety standards | EN 60601-1 |
| Electromagnetic compatibility | according to EN 60601-1-2 |

Operating conditions

| | |
|----------------------|-------------------------------------|
| Temperature range | of + 5 °C to + 40 °C |
| Humidity range | of 15% RH to 93% RH, non-condensing |
| Atmospheric pressure | 70 kPa to 106 kPa |

Storage and transport conditions

| | |
|----------------------|---------------------------|
| Temperature range | of MIN -25° C MAX +70° C |
| Humidity range | of MIN 10% RH –MAX 93% RH |
| Atmospheric pressure | 50 kPa to 106 kPa |

Replacement parts and optional accessories

| | |
|---|---------|
| Adult mask - large (ages 5 years plus) | 1127875 |
| Pediatric mask - medium (ages 1 – 5 years) | 1127798 |
| Pediatric mask - small (ages 0 – 18 months) | 1127822 |
| Mask adapter | 1125985 |
| Mouthpiece assembly | 1128501 |
| European (2 pin) plug adapter | 1127650 |
| UK (3 pin) plug adapter | 1127652 |
| Carry case | 1128576 |


Technical data EN13544-1

| | |
|--|---------------------------------|
| MMAD | 3.99 µm |
| Aerosol output | 1.19 mL |
| Aerosol output rate | 0.26 mL/min |
| Max fill volume | 8 ml |
| Max medication temperature increase at max fill | <10°C above ambient temperature |
| Noise level | <35dB at 1m |
| Distribution of particles >5 µm | 35.6% |
| Distribution of particles in range 2 - 5 µm | 52.3% |
| Distribution of particles in range <2 µm | 12.1% |
| Respirable fraction (%<5 µm) | 64.4% |
| Percentage of fill volume (2.5 mL) as aerosol output (delivered to filter) in one minute | 10.2% |
| Residual volume | 0.30 mL |

If the device has been stored at the extremes of the storage temperature, please allow at least 1 hour at room temperature before using the device.

Use of a face mask may affect the device performance. Data available on request.

Cumulative drug deposition – Salbutamol

Performance information provided as required by EN 13544-1:2007 may not apply for medications in suspension or high viscosity form. In such cases information should be sought from the medication supplier. Performance may vary based on atmospheric pressure depending on altitude above sea level, barometric pressure, and temperature.

Nebulizer performances are based upon testing that utilizes adult ventilatory patterns and are likely to be different from those stated for pediatric infant populations.

Electromagnetic information: The InnoSpire Go needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this document. Portable and mobile RF communications equipment can affect medical electrical equipment.

The CE mark on the product denotes compliance with all applicable EU Directives. Note that the Notified Body number does not apply to the RoHS (restriction of the use of certain hazardous substances in electrical and electronic equipment) Directive.

Expected service life: Nebulizer handset, battery and mask adapter, 3 years from date of purchase
 Mouthpiece assembly, 1 year from date of first use
 Masks, 1 year from date of first use

Symbol glossary

- | | | |
|---|--|--|
| <ul style="list-style-type: none"> ON (power) OFF (power) Serial number Follow instructions for use Complies with RTCA/DO-160F section 21, category M | <ul style="list-style-type: none"> CAUTION Alternating current Separate collection Atmospheric pressure Manufacturer Type BF applied parts | <ul style="list-style-type: none"> Class II internally powered device, double insulated Direct current Temperature limitation Humidity Date of manufacture |
|---|--|--|

Warranty

Respironics, Inc warrants the nebulizer handset and battery to be free from defects in materials and workmanship under normal use and operation for a period of 2 years from date of purchase from Respironics, Inc. The warranty is limited to repair or replacement, at Respironics, Inc's sole option, of any such component or equipment claimed to be defective when claim is shown to be bona fide by evaluation by Respironics, Inc. This warranty does not extend to any components or equipment subjected to misuse, improper operation, accidental damage, or unauthorized repairs, and is not extended to charges of, or for, labor repairs. All items returned must be properly packaged and shipped, prepaid, by the product distributor servicing the unit. Respironics, Inc shall not be liable to purchaser or others for loss of use of equipment or for indirect, or incidental or consequential damages that might arise.

Note: For Australian and New Zealand customers this warranty replaces the warranty contained above.

1. Respironics, Inc., a Philips Healthcare company ("Philips Respironics") warrants that the Products shall be free from defects of workmanship and materials and will perform in accordance with the Product specifications. 2. This warranty is valid for a period of 2 years from the date of purchase from an authorised Philips Respironics dealer. 3. If the Product is found to contain a defect of workmanship or materials or fails to perform in accordance with the Product specifications during the applicable warranty period, Philips Respironics will repair or replace, at its option, the defective material or part. 4. The customer is responsible for returning the product to an authorised Philips Respironics dealer, and collecting the product from the authorised Philips Respironics dealer after repair or replacement, at its own cost. Philips Respironics is responsible only for the freight cost of transporting the product between the authorised Philips Respironics dealer and Philips Respironics. Philips Respironics reserves the right to charge an evaluation and postage fee for any returned Product as to which no problem is found following investigation. 5. This warranty does not cover damage caused by accident, misuse, abuse, alteration, and other defects not related to materials or workmanship. 6. The warranty provided by Philips Respironics herein is not transferrable by the Buyer in the event of any sale or transfer of Products purchased by the Buyer from an authorised Philips Respironics dealer. 7. To the extent permitted by law, where the Buyer has the benefit of an implied guarantee under the Australian Consumer Law, but the Product is not of a kind ordinarily acquired for personal, domestic or household use or consumption Philips Respironics' liability shall be limited, at the option of Philips Respironics, to the replacement or repair of the Product or the supply of an equivalent Product. 8. To exercise your rights under this warranty, contact your local authorised Philips Respironics dealer. A list of all authorised dealers is available at the following link: <http://www.philips.com.au/healthcare>.

Alternatively, you can make a claim under this warranty by contacting Respironics directly at: Philips Electronics Australia Limited, 65 Epping Road, North Ryde NSW 2113, Australia. Tel: 1300 766 488, Email: prcontact@philips.com. 9. The following statement is provided to a Buyer who is a "consumer" under the Australian Consumer Law: Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the good repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. 10. The following statement is provided to a Buyer who is a "consumer" under the Consumer Guarantees Act 1993, New Zealand: Our goods come with guarantees that cannot be excluded under the Consumer Guarantees Act 1993. This guarantee applies in addition to the conditions and guarantees implied by that legislation.

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A. Instructions for use

1. The InnoSpire Go needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this document;
2. Portable and mobile RF communications equipment can affect the InnoSpire Go.

B. Technical description

1. Warning: the use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by the manufacturer of the InnoSpire Go as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the InnoSpire Go;
2. Warning: the InnoSpire Go should not be used adjacent to or stacked with other equipment.
- 3.

| Guidance and manufacturer's declaration – electromagnetic emissions | | |
|---|------------|---|
| The InnoSpire Go is intended for use in the electromagnetic environment specified below. The customer or the user of the InnoSpire Go should assure that it is used in such an environment. | | |
| Emissions test | Compliance | Electromagnetic environment – guidance |
| RF emissions CISPR 11 | Group 1 | The InnoSpire Go uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | The InnoSpire Go is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations/ flicker emissions IEC61000-3-3 | Complies | |

4.

| Guidance and manufacturer's declaration – electromagnetic immunity | | | |
|---|--|--|---|
| The InnoSpire Go is intended for use in the electromagnetic environment specified below. The customer or the user of the InnoSpire Go should assure that it is used in such an environment. | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment – guidance |
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±8 kV contact ±15 kV air | ±8 kV contact ±15 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient/burst IEC 61000-4-4 | ±2 kV for power supply lines | ±2 kV for power supply lines | The InnoSpire Go is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Surge IEC 61000-4-5 | ±1 kV line(s) to line(s) | ±1 kV line(s) to line(s) | |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11 | 0 % U _T (100 % dip in U _T) for 0.5 cycle 0 % U _T (100 % dip in U _T) for 1 cycle 70 % U _T (30% dip in U _T) for 25 cycles (50Hz) for 30 cycles (60Hz) 0 % U _T (100% dip in U _T) for 250 cycles (50Hz) for 300 cycles (60Hz) | 0 % U _T (100 % dip in U _T) for 0.5 cycle 0 % U _T (100 % dip in U _T) for 1 cycle 70 % U _T (30% dip in U _T) for 25 cycles (50Hz) for 30 cycles (60Hz) 0 % U _T (100% dip in U _T) for 250 cycles (50Hz) for 300 cycles (60Hz) | |
| Power frequency (50/60Hz) magnetic field IEC 61000-4-8 | 30 A/m 50Hz or 60Hz | 30 A/m 50Hz or 60Hz | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |
| NOTE U _T is the a.c. mains voltage prior to application of the test level. | | | |

5.

| Guidance and manufacturer's declaration – electromagnetic immunity | | | |
|---|---|------------------|---|
| The InnoSpire Go is intended for use in the electromagnetic environment specified below. The customer or the user of the InnoSpire Go should assure that it is used in such an environment. | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment – guidance |
| Conducted RF IEC 61000-4-6 | 3Vrms 150kHz to 80MHz outside ISM bands ^a | 3Vrms | Portable and mobile RF communications equipment should be used no closer to any part of the InnoSpire Go, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance d = 1.2√P d = 2.0√P |
| | 6Vrms 150kHz to 80MHz in ISM and amateur radio bands ^a | 6Vrms | |
| Radiated RF IEC 61000-4-3 | 10V/m 80MHz to 2.7GHz | 10V/m | d=1.2√P 80MHz to 800MHz d=2.3√P 800MHz to 2.7GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m) ^b . Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^c , should be less than the compliance level in each frequency range ^d . Interference may occur in the vicinity of equipment marked with the following symbol:  |
| <p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> <p>a. The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.</p> <p>b. The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.</p> <p>c. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the InnoSpire Go is used exceeds the applicable RF compliance level above, the InnoSpire Go should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the InnoSpire Go.</p> <p>d. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p> | | | |

6.

| Recommended separation distances between portable and mobile RF communications equipment and the InnoSpire Go | | | | |
|--|--|--|------------------------------|-------------------------------|
| The InnoSpire Go is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the InnoSpire Go can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the InnoSpire Go as recommended below, according to the maximum output power of the communications equipment. | | | | |
| Rated maximum output power of transmitter W | Separation distance according to frequency of transmitter m | | | |
| | 150 kHz to 80 MHz outside ISM bands d=1.2√P | 150 kHz to 80 MHz inside ISM bands d=2.0√P | 80 MHz to 800 MHz d=1.2√P | 800 MHz to 2.7 GHz d=2.3√P |
| 0.01 | 0.12 | 0.20 | 0.12 | 0.23 |
| 0.1 | 0.38 | 0.63 | 0.38 | 0.73 |
| 1 | 1.2 | 2.0 | 1.2 | 2.3 |
| 10 | 3.8 | 6.3 | 3.8 | 7.3 |
| 100 | 12 | 20 | 12 | 23 |
| For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. | | | | |
| NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. | | | | |
| NOTE 2 The ISM (Industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. | | | | |
| NOTE 3 An additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.7 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. | | | | |
| NOTE 4 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. | | | | |

© 2017 Koninklijke Philips N.V. All rights are reserved.
Trademarks are the property of Koninklijke Philips N.V. or their respective owners.
Respironics Inc. and Respironics Respiratory Drug Delivery (UK) Ltd are Philips companies.

CE0123

Manufactured for
Respironics Respiratory Drug Delivery (UK) Ltd
Chichester Business Park, City Fields Way
Tangmere, Chichester PO20 2FT UK

www.philips.com/respironics

PN 1126300 R03