




# ABPMpro

**A multi-sensor recording unit, including the measurement of ambulatory (24h) blood pressure**



## Technical design

- Oscillometric sensor
- 1-channel-ECG sensor (cuff)
- Actigraphy sensor (X/Y/Z)
- Body position sensor
- PPG (optical) sensor (ext)
- 3-channel-ECG sensor (ext)
  
- 125 gram weight (4.4 oz)  
incl. battery, pump, valve, etc.




## Device configurations

### ABPMpro Classic:

- Oscillometric (cuff) blood pressure
- Accelerometry (tri-axial)
- Body position
- ECG heart rate variability (internal)


ABPMpro Classic



### ABPMpro Cardio (adding 3-lead ECG cable):

- 3-channel ECG recording (external)
- Respiration rate
- Pre-ejection period


ABPMpro Cardio



### ABPMpro Research (adding ECG and PPG):

- PTT beat-to-beat blood pressure
- Sleep arousals

ABPMpro Research



## Device configurations

### ABPMpro Classic:


Oscillometric (cuff) blood pressure

Accelerometry (tri-axial)

Body position


ECG heart rate variability (internal)

ABPMpro Classic




## Oscillometric (cuff) blood pressure

Measurement during **cuff inflation** and/or  
Measurement during **cuff deflation**



**artefact during inflation**

**2nd attempt during deflation**



## Validation data

The oscillometric measurements of blood pressure has been validated according to the **ISO 81060-2** standard for both, **inflation and deflation** recordings

Including the validation post **change in heart rate** within the same individuals (ambulatory condition)

## Validation of the ABPMpro ambulatory blood pressure monitor in the general population according to AAMI/ESH/ISO Universal Standard (ISO 81060-2:2018)

Bernhard Roth<sup>a</sup>, Tomas Lucca Bothe<sup>b</sup>, Andreas Patzak<sup>b</sup> and Niklas Pilz<sup>b</sup>

**Conclusion** The ABPMpro device fulfilled the ISO 81060-2:2018 requirements in the general population and in the ambulatory setting and can therefore be recommended for clinical use. *Blood Press Monit XXX:*

Correspondence to Andreas Patzak, MD, PhD, Institute of Translational Physiology, Charite – Universitätsmedizin Berlin, corporate member of Freie Universität Berlin, Humboldt-Universität zu Berlin, and Berlin Institute of Health, Chariteplatz 1, 10117 Berlin, Germany  
Tel: +49 30 450 528220; fax: + 49 (0)30 4507528972;  
e-mail: andreas.patzak@charite.de

## Validation data



JOINT INITIATIVE WITH



---

## VALIDATED DEVICES FOR AMBULATORY BLOOD PRESSURE MONITORING

### Preferred devices (upper-arm) (21)

---

**A&D** TM-2440 \*

**A&D** TM-2441

**Andon** iHealth CardioMed ABP100

**Beneware** ABP-021

**Custo Med** custo screen 400

**Custo Med** custo screen pediatric

**EDAN** SA-10

**Hingmed** WBP-02A

**KANG** KC-2850

**Meditech** ABPM-06

**Microlife** Watch BP O3

**Microlife** WatchBP O3 (BP3SZ1-1) \*

**Microlife** WatchBP O3 AFIB \*

**Novacor** Diasys 3 (DIS-0001-00)

**Novacor** Diasys 3 Plus (DIP-0001-00)

**PAR Medizintechnik & Co.** PHYSIO-PORT UP

**PAR Medizintechnik & Co.** TONOPORT VI

**Philips** DL8760

**Somnomedics** ABPMpro

**SpaceLabs** 90227


**TaiDoc** TD 3127AT Uright

## Accelerometry (tri-axial)

- Allows to trigger **an acoustic signal** and/or to **slightly delay** the scheduled cuff measurement, when the sensor records too much movement (at the time of scheduled recording)
- Simultaneous recording of blood pressure and patient's activity data supports the **interpretation of blood pressure changes**

peak in blood pressure


increased motoric activity





## Body position

Body position sensor together with the data on motoric activity allow for a more **precise identification of day-/ versus night-time period** (in-bed time)




**Recorded body positions: supine, prone, upright, left, right, walking**

**-> Objective data, not dependent on the memory of the patient.**

## ECG heart rate variability (internal)

**Cuff-integrated ECG electrodes** allow for the continuous recording of 1-channel-ECG data and subsequent analysis of heart rate variability



## Device configurations


### ABPMpro Classic:

- Oscillometric (cuff) blood pressure
- Accelerometry (tri-axial)
- Body position
- ECG heart rate variability (internal)


### ABPMpro Cardio (adding 3-lead ECG cable):

- 3-channel ECG recording (external)
- Respiration rate
- Pre-ejection period

ABPMpro Classic




ABPMpro Cardio





## 3-channel ECG recording (external) Respiration rate Pre-ejection period

Attaching a **3-lead ECG cable** to the ABPMpro monitor allows for the simultaneous recording of 24h holter ECG data.

With the use of **impedance cardiography**, respiration rate and pre-ejection period are derived



ABPMpro Cardio




## Device configurations

### ABPMpro Classic:

Oscillometric (cuff) blood pressure  
 Accelerometry (tri-axial)  
 Body position  
 ECG heart rate variability (internal)


ABPMpro Classic



### ABPMpro Cardio (adding 3-lead ECG cable):

3-channel ECG recording (external)  
 Respiration rate  
 Pre-ejection period


ABPMpro Cardio



### ABPMpro Research (adding PPG):

PTT beat-to-beat blood pressure  
 Sleep arousals

ABPMpro Research




## Sleep arousals

Attaching a **photoplethysmogram (PPG)** sensor to the ABPMpro monitor allows for the simultaneous recording of optical pulse waves.

The reduction in **pulse wave amplitude** is an accepted method to detect autonomous arousals

**Detecting arousals by use of the PPG signal on the upper-arm is an indicator for sleep fragmentation**





Age	Normvalue
Young and middle aged adults	15 / h
> 60 Years	27 / h


## PTT beat-to-beat blood pressure

Attaching a **photoplethysmogram (PPG)** sensor to the ABPMpro monitor allows for the simultaneous recording of optical pulse waves.

Pulse transit time (PTT) = Pulse traveling time from **proximal** (ECG r-peak) to **distal** (optical pulse wave) site, corrected for pre-ejection-period




**Changes in PTT translate into changes of blood pressure. An increase in PTT reflects a decrease in blood pressure and vice versa.**



## PTT beat-to-beat blood pressure

- Continuous blood pressure (beat-to-beat) based on PTT/PAT (pulse transit time / pulse arrival time)

**Oscillometric recording at 9:05 pm**  
106 / 66 mmHg



**Oscillometric recording at 9:20 pm**  
121 / 78 mmHg

**PTT/PAT recording 15 minute beat-to-beat BP recording**

**Each cuff-BP (24h) used for re-calibration of PTT/PAT**




## Limitations

The cuff measurement of ABPMpro is obtained via the **oscillometric** measurement principle. As such, the same limitations as for all automated cuff blood pressures apply, including:

- **Movement artefacts**
- **Severe arrhythmia**
- **Left-ventricular assist devices ('heart pumps')**
- **Measurement during hemodialysis (blood volume changes)**
  
- **BP ranges (systolic) below 60 mmHg / above 230 mmHg**
- **BP ranges (diastolic) below 40 mmHg / above 130 mmHg**
  
- **Validation in children (< age 12) pending**
- **Validation in pregnancy pending**
- **Validation in atrial fibrillation pending**
- **Validation of beat-to-beat blood pressure pending**

## Plotting the curve – opening the black box



Plotting the pressure curves is a helpful tool to **eliminate cuff errors**



## Plotting the curve – opening the black box



Tri-axial accelerometer to **eliminate movement artefacts** and identifying the magnitude of influence of movement on the pressure curve





## Technical features (2): Measurements

- Time-stamped multi-sensor data
  - Many more channels to chose from for graphical display
- Oscillometric (cuff) blood pressure
  - PTT beat-to-beat blood pressure
  - 3-channel ECG recording
  - Actigraphy / accelerometry / gyroscopy
  - Body poistion (arm/thorax)
  - Sleep arousals
  - Respiration rate
  - Heart rate variability
  - Pre-ejection period



## Technical features (3): Measurements

- Time-stamped multi-sensor data  
(track what happens)

**Irrgeularity in the ECG, ...**

**... will also show in the PPG**



## Technical features (4): Export all data

### High sampling rate, ...

Patient Information Data

Identification    Diagnosis    **Montage**    Logbook    Track

Device ID:     Firmware Version:     Cuff size:

Day: 6 h (15 min.), Night: 22 h (30 min.), Maximum: 180 mmHg, Autostop

Abort Recording:

Number	Description	Signal Type (Port)	Rec. Rate (Hz)
0	Accu	Accu	4
1	Pt. Marker	Patientenmarker	4
2	X	X	32
3	Y	Y	32
4	Z	Z	32
5	Cuff press	Cuff press.	128
6	ECG	ECG	256
7	ECG I	ECG I	256
8	ECG II	ECG II	256
9	ECG III	ECG III	256
11	Imp-Resp	Imp-Resp	256
12	Position	Pos. Ext	4
13	Pleth raw	Pleth Raw	256
14	Log	Log	256

OK    Cancel

## ... export of all raw data!

Export

Export to own subfolder

Export Path:  ...

Filename:

RIFF Export     ASCII Export

Complete Recording     TIB area

Selection     User Areas

without Time Base

Channel name	Export rate
<input checked="" type="checkbox"/> ECG II	256 Hz
<input checked="" type="checkbox"/> Pleth	256 Hz
<input checked="" type="checkbox"/> Cuff press	128 Hz
<input checked="" type="checkbox"/> Cuff osc.	128 Hz
<input checked="" type="checkbox"/> Activity	32 Hz
<input checked="" type="checkbox"/> ECG	256 Hz
<input checked="" type="checkbox"/> ECG I	256 Hz
<input checked="" type="checkbox"/> X	32 Hz
<input checked="" type="checkbox"/> Y	32 Hz
<input checked="" type="checkbox"/> Z	32 Hz
<input checked="" type="checkbox"/> Accu	4 Hz
<input checked="" type="checkbox"/> ECG III	256 Hz
<input checked="" type="checkbox"/> Imp-Resp	256 Hz
<input checked="" type="checkbox"/> Position	4 Hz

Export    Cancel



## Technical features (5a):

Export all data – data format in ASCII / RIFF / EDF

### Activity

Signal Type: Mag\_Type  
 Start Time: 28.02.2022 12:19:41  
 Sample Rate: 32  
 Length: 2764800  
 Unit: mg

Data:  
 508  
 118  
 121  
 131  
 129  
 133  
 132  
 132  
 128  
 131  
 140  
 141  
 138  
 135  
 133  
 ...  
 ...  
 ...

### Cuff osc.

Signal Type: BD\_AC\_Type  
 Start Time: 28.02.2022 12:19:41  
 Sample Rate: 128  
 Length: 11059200  
 Unit: mmHg

Data:  
 0,0  
 0,0  
 0,0  
 0,0  
 0,0  
 0,0  
 0,0  
 ...  
 ...  
 ...  
 0,01  
 0,07  
 0,13  
 0,2  
 0,27  
 0,34  
 ...  
 0,42  
 0,49  
 0,55

### Cuff press.

Signal Type: BD\_Type  
 Start Time: 28.02.2022 12:19:41  
 Sample Rate: 128  
 Length: 11059200  
 Unit: mmHg

Data:  
 0,0  
 0,0  
 0,0  
 0,0  
 0,0  
 0,0  
 0,0  
 0,0  
 0,0  
 ...  
 ...  
 ...  
 127,45  
 127,38  
 127,37  
 127,37  
 127,33  
 127,3  
 127,32

### ECG channel #1

Signal Type: ECG1\_Type  
 Start Time: 28.02.2022 12:19:41  
 Sample Rate: 256  
 Length: 22118400  
 Unit: µV

Data:  
 986  
 958  
 997  
 1009  
 1074  
 1143  
 1021  
 927  
 856  
 886  
 ...  
 ...  
 ...  
 25  
 25  
 22  
 20  
 14





## Technical features (5b):

Export all data – data format in ASCII / RIFF / EDF

### Imp-Resp

Signal Type: ECG\_Imp\_Resp\_Type  
Start Time: 28.02.2022 12:19:41  
Sample Rate: 256  
Length: 22118400  
Unit: Ohm

Data:

1132,4  
1132,5  
1132,4  
1132,4  
1132,4  
1132,4  
1132,4  
1132,5  
1132,5  
1132,5  
...  
...  
...  
1115,1  
1115,1  
1115,1  
1115,1  
1115,1  
1115,0  
1115,1  
1115,1  
1115,1

### Pleth

Signal Type: Pleth\_Type  
Start Time: 28.02.2022 12:19:41  
Sample Rate: 256  
Length: 22118400  
Unit:

Data:

...  
...  
...  
2796  
2802  
2808  
2814  
2820  
2826  
2831  
2836  
2841  
2846  
2851  
2855  
2859  
2863  
2866  
2869  
2872

### Position

Signal Type: Pos\_XYZ\_Type  
Start Time: 28.02.2022 12:19:41  
Sample Rate: 4  
Length: 345600  
Unit:

Data:

Upright  
Upright  
Upright  
Upright  
...  
...  
...  
Prone  
Upright  
Upright  
Upright  
Upright  
Prone  
Prone  
Prone  
Prone  
...  
...  
...  
Supine  
Supine


### X-axis

Signal Type: KLX\_Type  
Start Time: 28.02.2022 12:19:41  
Sample Rate: 32  
Length: 2764800  
Unit: mg

Data:

-868  
-867  
-866  
-866  
-866  
-866  
-865  
...  
...  
154  
155  
155  
154  
154  
...  
...  
-708  
-778  
-916  
-1000

**In addition (not shown): ECG channels #2 and #3, battery voltage, Y- and Z-axis**



**ABPMpro**

## Comparison to standard ABPM devices (page #1)

ABPMpro Classic

Characteristic	ABPMpro	Standard ABPM
<b>No tube design</b> (device attached to cuff)	✓	✗ (✓)
<b>Inflation-/deflation measurement</b>	✓	✗
<b>Actigraphy / accelerometry sensors</b>	✓	✗ (✓)
<b>Body position</b> (arm sensor)	✓	✗ (✓)
<b>Internal ECG</b> inside the cuff for heart-rate-variability	✓	✗
<b>Validated</b> ISO81060:2	✓	✗ (✓)

## Comparison to standard ABPM devices (page #2)

Characteristic	ABPMpro	Standard ABPM
<b>Continuous blood pressure</b> (beat-to-beat)	✓	✗
<b>Optical sensor</b> for detection of sleep arousals	✓	✗
<b>3- channel ECG</b> recording (holter ECG)	✓	✗ (✓)
<b>Body position</b> (thorax sensor)	✓	✗ (✓)
<b>Impedance cardiography</b> for detection of respiration rate	✓	✗
<b>Export of all raw data</b> incl. synchronized time stamp	✓	✗

## Take home messages

Improved **patient acceptance** to ABPM recordings  
(increased patient comfort)



More successful (valid) ABPM recordings over 24h – **less repeat exams**  
(avoiding tube kinks and „two-shots-on-one-recording“)

Additional data for the **interpretation of 24h blood pressure** profiles  
(activity data, precise time-in-bed calculation)

ABPMpro research: Simultaneous recording of BP, ECG, activity, PPG, etc.  
and the interaction among these parameters may support in the  
understanding of underlying **mechanisms of intervention**

**Export of all raw signal data** will allow for data post-processing and the  
application of machine learning and artificial intelligence