**HERSILL**

Medical Devices

**Medical Suction Equipment**

*Directions for Use and Maintenance*

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 **HERSILL**  
www.hersill.com

## 1. Precautions and Warnings

**CAUTION:** Federal law (USA) restricts this device to sale by or on the order of a physician (CFR 21 801.109).

- The V7 unit must be used by trained personnel in the use of medical equipment.
- Read and follow the instructions of this guide.

**Risk of explosion or fire:** Do not use the machine near of flammable or explosive liquids or gases.

**Risk of electric shock:** Do not submerge this machine under any liquid. It could cause an electric discharge hurting to the people or damage the machine itself. Unplug the equipment before cleaning the unit. Use a minimum quantity of water or other liquids to avoid electric discharge.

Do not remove cover. Refer servicing to qualified personnel.

To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth

**Risk in extreme environmental conditions:** Do not use this machine under environmental conditions that exceed the specified limits: 0-55°C and 95% humidity; this could be dangerous to the security and affect the functioning in a negative way.

**Risk of liquid entrance inside the pump:** Do not start the machine if there is any deterioration in the overfilling system of the collector jar. Do not start the machine if there exists or if you suspect the existence of contamination owing to excess liquid. Return the device for service

**Risk of contaminating the patient:** Always use a clean collector jar and a new bacteria filter for each patient. After each use, it is necessary to disinfect the whole collector jar (see *Chapter 7.2*)

**Risk for the machine owing to a wrong connection:** While changing the battery, make sure not to exchange the polarity, the poles have the same size and it is possible to connect them incorrectly (see *Chapter 7.1.2*).

- Remove the battery when the device is not used for long periods of time.
- To insure proper grounding reliability, connect to "Hospital grade" receptacle or equivalent.
- Repairs are to be carried out by personnel qualified by HERSILL or one of its authorized dealers. This device must be clean of all blood or other organic materials prior to returning for service.
- Do not modify this equipment without authorization of the manufacturer.
- Please follow the instructions on the use of cleaning products and disinfectants.
- Please use only the accessories provided by HERSILL or its authorized distributors.
- Old electrical and electronic equipment must be disposed separately and may not be included in regular domestic waste. Disposal of device, components and accessories must strictly conform to local laws and regulations as set out by the relevant authorities.
- The configuration of the components of the product may be subject to variations without previous notification. We recommend that you consult with the HERSILL dealer to find out the latest version.

## 2. Introduction

### 2.1. Indications for use

The V7 unit is a electrically powered portable suction equipment, intended for use in clinical institutions; the V7 m, V7 mx, V7 plus b, V7 plus b emergency models are specifically designed for transport and emergency field use.

It is indicated for taking out secretions, blood or vomits from the respiratory track of the patients to let air flow freely to the lungs.

There are a full range of V7 portable aspirators. There is a specific machine that adapts to each need as far as the vacuum and flow, the collection capacity, different modes and levels of electrical supply are concerned.

### 2.2. Contraindications

None known

### 2.3. Conformity with product standards and regulations

ISO 10079-1:1999

EN IEC 60601-1:2006

EN IEC 60601-1-2:2007

EN IEC 60601-1-6:2010

EN 1789:2007+A1:2010

## 3. Symbols

<b>IPX1</b>	Degree of protection provided by enclosures against dripping water
	Class II (IEC 60601-1)
	Type BF (IEC 60601-1)
	Direct current (DC)
	Attention. Consult accompanying documents
	CE Mark. This product is in compliance with MDD 93/42/EEC

## 4. Description of the equipment

### 4.1. General description

- A. Main unit (casing)
- B. Liquids collector bottle
- C. Suction tubes (must be changed for each patient).
- D. Bacteria filter (must be changed for each patient).
- E. Electrical connection cables (depending of the models)
- F. ON/OFF switch
- G. Vacuum indicator (vacuum-meter)
- H. Vacuum regulation control knob
- I. LED indicators (only for model *V7 m*, *V7 mx*, *V7 plus b* and *V7 plus b emergency*)

---

### 4.2. Led indicators and vacuum control knob

Power on LED indicator

Vacuum control knob

External power LED indicator

Battery status LED indicators

## ENGLISH

Power ON LED indicator:

This green LED is continuously lit while the V7 unit is switched ON.

External power LED indicator:

This green LED is continuously lit while external AC- or DC-power is connected.

Battery status LED indicators:

3 LEDs: fully charged, 75-100% battery capacity



2 LEDs: 50-75% battery capacity

1 LED: 25-50% battery capacity

 Device settings, state of battery, temperature, etc can influence battery status

Vacuum control knob: Turn right: increase vacuum  
Turn left: decrease vacuum

### 4.3. Electrical connection cables

## 5. Preparing of the equipment for its use

The V7 suction unit is supplied assembled with the exception of the lowest suction tubes to allocate the bacteria filter between the suction unit and the collection bottle; afterwards the equipment will be ready to be used. (see Chapter 5.1.).

The collector jar, the suctioning tubes and the bacteria filter are the only components that have to be changed for each patient and are easily removable. On assembling the equipment, special care has to be taken to respect the direction of the suctioning flow that is identified and recorded in the collector jar.

Never operate the unit without bacterial filter, the unit could be damaged

Do not forget to check the correct operation voltage

### 5.1. Previous assembly procedure

1. Connect the 5cm tube between the bacteria filter and the input of the suction unit.
2. Connect the 15cm tube between the bacteria filter and the outlet *VACUUM* of the collection bottle.
3. Put the suction tube (1.5mx8mm) between the outlet *PATIENT* of the collection bottle and the suction terminal and roll it in the tube location of the unit.
4. Connect the electric power cable as required by the model.
5. Switch on the unit (ON/OFF switch).
6. Collapse the flow twisting the suction tube and checking the vacuum meter to test that the unit makes vacuum without any leakage and that this level varies when the vacuum regulator is acted on.
7. (Only in models with battery) Check the battery level while in function; in case of need, leave it connected to a power supply during the time needed (the full charge of a new battery that is run down is five hours).
8. Switch off the unit (ON/OFF switch).

### 5.2. Final assembling for a specific patient

1. Connect the suction terminal to a suction catheter specific for this application and patient.
2. Switch-on the unit (ON/OFF SWITCH).
3. Obstruct the flow in the suction catheter and adjust the level of vacuum prescribed in the vacuum-meter turning the vacuum regulator control knob.
4. The unit is ready to put the suctioning catheter to the patient.

## 6. Functioning and Operation

-  Do not start the machine if you suspect the presence of explosive or inflammable liquids or gases, water inside the machine or environmental conditions that are not within the specified limits.
-  Do not start the unit if the liquid level of the jar exceeds 50% of its capacity.
-  Do not start the machine when the battery charger is connected to the mains. (V7 Plus b)
-  When operating unit in battery mode, the user must monitor the LED's to ensure enough battery power is left to finish procedure.

- Verify the assemblage and test it as detailed above in *Chapter 5*.
- Proceed to a visual inspection: check that there is no one damaged component and the unit and the collection jar are clean.
- Put the unit in a straight and flat area.
- Apply the suction catheter to the patient.
- When the liquid level reaches 50% of the total volume of the collection bottle, it is recommended to empty or change it before the safety valve acts, in this way the risk of the entrance of humidity in the bacteria filter is minimized.
- Once used, disassemble suction line and proceed to its clearing / substitution as indicated in *Chapter 7*. After this will be necessary to assemble the unit as detailed in *Chapter 5*.

## 7. Maintenance and Cleaning

### 7.1. Periodic (planned) maintenance

Clean the unit after every use or maintenance.

Perform an operation test before each use and after any maintenance operation and at least once a month (see *Chapter 5.2*).

In the models *V7 m*, *V7 mx*, *V7 plus b* and *V7 plus b emergency*, do a battery capacity test every 6 months (see *Chapter 7.1.1*) and charge the battery for 5 hours after each use without connecting it to the mains, and every 2 weeks.

It is recommended to replace the battery every 3 years, or before in the case to detect any fail in the scheduled tests (see *Chapter 7.1.2*).

#### 7.1.1. Test of capacity of the battery

##### **(V7m, V7mx, V7plusb and V7plusb emergency)**

1. Charge the battery for 6 hours and disconnect from the supply.
2. Switch-on the unit and turn the vacuum regulator to high vacuum, make sure that the vacuum indicator is at the correct position.
3. Suction speed: Measure the time needed to suck 500 ml of water.
4. Repeat this step 5 times; the average time should be less than 10 seconds.
5. Let the unit work continuously with free flow for 20 minutes.
6. If the unit stops before 20 minutes or it does not pass the suction speed test, substitute the battery with a new one pursuant Chapter 7.1.2.

### 7.1.2. Battery replacement

*(V7m, V7mx, V7plusb and V7plusb emergency)*

1. Unplug the unit from the mains.
2. Withdraw the 6 screws from the back side.
3. Maintain the unit in vertical position and separate the both parts of the casing.
4. Take off the old battery terminals and recycle it pursuant the local rules.
5. Set up the new battery with its terminals in its correct polarity.
6. Close the casing and set up the 6 screws in the back-side.
7. Charge the battery for 5 hours and test its capacity (*Chapter 7.1.1*).

### 7.1.3. Fuse replacement

Unplug the unit from the mains

Replace the fuse with a new one of the same rating. (See Technical specifications)

12/24Vdc fuse:  
turn left using a  
flat head  
screwdriver and  
take the fuse  
holder out

110-240Vac fuse:  
Using a flat head  
screwdriver, pry out  
the plastic fuse holder

## 7.2. Cleaning and Disinfection

*Disconnect the machine from any external mains before cleaning it. Use a minimum quantity of water or other liquids to avoid an electric discharge.*

*Contaminated devices must be thoroughly cleaned prior to disinfection or sterilization*

*Personal protective equipment must be worn when handling contaminated equipment, particularly tubing, filter and collection jar (protective clothing, disposable latex gloves, mask, eye protection)*

*If any liquid accidentally enters in the main unit or in the pump the unit may be damaged. For this reason, we recommend not to use the aspirator itself as a method to clean the collection jar (taking in water with disinfectants) nor connecting the suctioning tubes in the opposite way (the safety device would not work). The failure should invalidate the guarantee and the technical specifications under the manufacturer's responsibility.*

*In the event that contaminated fluids have moved past the internal unit, remove the tubing and the filter, clean the casing (7.2.1), and decontaminate the collection bottles (7.2.2). Return the device for service.*

*Polycarbonate jars can be reused up to 40 times and Polysulfone jars up to 100 times. Before each use, check that collection jars are not damaged.*

*To obtain additional information including questions on infection control procedures, please contact Hersill or your dealer. Telephone: +34 91 616 41 11*

### 7.2.1. Cleaning of the main unit (casing)

1. Use a cloth or sponge that is dampened with a mild detergent to clean the external surfaces of the main unit.
2. Repeat the previous step with a mild disinfectant compatible with ABS and PC.
3. Use a cloth or sponge dampened with water and wipe the surfaces again. Dry the surfaces using a clean cloth or a paper towel.

⚠ Never use dissolvent or abrasive products.

### 7.2.2. Decontamination of the liquids collection bottle

- ⚠ It is the user's responsibility to qualify any deviations from the recommended method of processing.
- ⚠ Follow the decontamination product's instructions. Check the decontamination product compatibility with the materials
- ⚠ After decontamination, label products with the date
- ⚠ The disposable jars should be processed according to the local rules.


After each use, decontaminate the collector recipient pursuant the next procedure:

**A. Disconnect** the 3 tubes and the filter and replace them for new ones.

**B. Disassemble** the cover and the floating device from the transparent jar:

Free the two locking crabs to open the jar, lift the cover and pull down the floating device.

If dirt is observed in the rubber quad-rings of the floating device or in the locking crabs, they can be disassembled (it is recommended to warm up in water at 50°C to expand the rubber rings).

**C. Washing and rinsing (manual or mechanical):**

Wet in lukewarm water (<43°C) before adding any detergent; Wash parts thoroughly in warm water using a compatible detergent (rub with brushes, if necessary); Rinse thoroughly with sterile water, by immersing it completely in three separate copious volumes of water. Each rinse should be a minimum of one minute. Do not reuse the water. Finally, allow to dry.

- Alkaline or neutral detergents: non-ionic tension-actives concentration lower than 15%, or ionic tension-actives concentration lower than 30%.
- Enzymatic detergents: enzymes concentration lower than 15%.

**D. Disinfection / Sterilization:**

Disinfected/Sterilized jars should be dated and stored in clean and dry areas to avoid recontamination

Avoid to put in contact the parts with the autoclave chamber, because the material can be damaged. It is recommended to use a basket.

All parts can be processed pursuant one of the next 2 methods:

- **Steam Sterilization in gravity autoclave:**  
 Polycarbonate (PC) jars: (max. 40 times)  
     at 121°C and 1.3 kg/cm<sup>2</sup> for 25 minutes, drying time > 30 minutes.  
 Polysulfone (PSU) jars: (max.100 times)  
     at 132-136°C and 2.4 kg/cm<sup>2</sup> for 15 minutes, drying time > 30 minutes.
- **High level Disinfection:**  
 With Glutaraldehyde (i.e. *Cidex*<sup>®</sup>): 2% solution at ambient temperature during 60 minutes. Rinse thoroughly with sterile water, by immersing it completely in three separate copious volumes of water. Each rinse should be a minimum of one minute. Do not reuse the water. Finally, allow to dry.

**E. Inspection, reassembly and perform the function test pursuant Chapter 5.**

**8. Troubleshooting information**

FAULT	CONDITION	ACTION
The V7 unit does not operate	Control Unit not turned on	Turn Main Switch „ON/OFF”
	Wrong operating voltage	Check mains voltage output
	Mains not connected	Connect mains to the Unit
	External power indicator lit	Send the unit to repair
	Defective Fuse	Check fuse and replace it
The V7 unit does not work from its battery	Battery status LEDs is not lit	Place the V7 unit on charge. If does not operate after charging, replace the battery
The V7 unit operates, but with little flow or little vacuum	Floater valve is closed	Unplug the Filter/Vacuum tubing to release the vacuum
	Collection jar full	Remove and replace the jar
	Defective connection from the Unit to the jar	Install the Vacuum Connector Tubing correctly
	Patient Suction Tubing twisted or blocked	Untwist, clear or replace the tubing

**9. EMC Information**

V7 equipment meets the IEC 60601-1-2 EMC requirements for medical devices.

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the accompanying documents. Portable and mobile RF communications equipment can affect medical electrical equipment.

## 10. Technical Specifications

	V7 m	V7 mx	V7 plus b emerg.	V7 plus b	V7 plus dc	V7 dc	V7 plus ac	V7 ac
<b>RANGE</b>	High Vacuum - High Flow	High Vacuum - High Flow	High Vacuum - High Flow	High Vacuum - High Flow	High Vacuum - High Flow	High Vacuum - Low Flow	High Vacuum - High Flow	High Vacuum - Low Flow
<b>SUPPLY</b>	110 - 220 Vac - 12/24 Vdc - Battery	100- 240 Vac - 12 Vdc - Battery	12 Vdc - Battery	12 Vdc - Battery	12 Vdc - Battery	12 Vdc - Battery	110 Vac 60 Hz or 220 Vac 50 Hz	110 Vac 60 Hz or 220 Vac 50 Hz
<b>VACUUM</b>	84 kPa 630 mmHg	84 kPa 630mmHg	84 kPa 630 mmHg	84 kPa 630 mmHg	84 kPa 630 mmHg	75 kPa 562mmHg	84 kPa 630 mmHg	80 kPa 600 mmHg
<b>FLOW</b>	30 L/min	30 L/min	30 L/min	30 L/min	30 L/min	14L/min	30 L/min	20L/min
<b>POWER</b>	130 W(AC) 84 W (DC)	160 W (AC) 84 W (DC)	85 W	85 W	85 W	18 W	120 W	250 W
<b>BATTERY</b>	12 V 4,5 Ah. Ac-Pb	12 V 4,5 Ah. Ac-Pb	12 V 4,5 Ah. Ac-Pb	12 V 4,5 Ah. Ac-Pb				
<b>AUTONOMY (battery)</b>	65 min	65 min	65 min	65 min				
<b>BATTERY CHARGER</b>	Internal (integrated in the unit)	Internal (integrated in the unit)	External (optional accessory)	External (optional accessory)				
<b>10 g tested wall bracket</b>			YES					
<b>WEIGHT</b>	5,2 kg	4,5 kg	4,2 kg	4,2 kg	2,5 kg	2,2 kg	3,5 kg	3 kg
<b>Fuse</b>	AC F2AL250V - DC F8AL250V	AC F2AL250V - DC F8AL250V	DC F8AL250V	DC F8AL250V	DC F8AL250V	DC F8AL250V	AC F5AL250V	AC F5AL250V
<b>Dimensions (packaged)</b>	460 x 210 x 290 mm							
<b>Noise level</b>	49 ± 1,5dB (AS) / 1 m							
<b>References (with jar 1L PC)</b>	5320231	5320236	5320281	5320261	5320251	5320211	5320241 (220 V) 5321241 (110V)	5320201 (220 V) 5321201 (110 V)

## 11. Accessories and Spare parts

<i>Description</i>	<i>Reference</i>
1L PC jar , with cap	4383001
1L PSU jar , with cap	4383000
1,7L PC jar , with cap	4383002
1,7L PSU jar , with cap	4383012
10 g wall bracket with 12 Vdc supply	4320300
10 g wall bracket (without 12 Vdc supply)	4320301
Basic wall bracket	4320290
Battery charger 2.5 A / IEC 60320/C8	4320033
Power cable 12/24Vdc	4320032
Battery 4,2A / 12V	3320028
Bacteria filter	3320534
PVC tube 1m 8x13mm	0940023
Transport cover-bag	3320223
Bracket for jars "Abbott"	4383015
Bracket for disposable jars	Consult

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<b>CE DECLARATION OF CONFORMITY</b>	
HERSILL, S.L. Puerto de Navacerrada, 3 - P.I. las Nieves 28935 Móstoles (Madrid) Spain	
declare under our sole responsibility that the products listed below are in conformity with the provisions of Directive 93/42/EEC of 14 June 1993 and its revised version, concerning medical devices.  <b><i>This declaration of conformity is subject to the procedure defined in the Annex II of the European Directive 93/42/EEC and its revised version.</i></b>	
<i>Product type</i>	<b>ELECTRICALLY POWERED MEDICAL SUCTION EQUIPMENT</b>
<i>Product name</i>	<b>V7 m , V7 mx , V7 plus b , V7 plus b emergency , V7 dc , V7 plus dc , V7 ac , V7 plus ac</b>
<i>Product Classification</i>	<b>II a</b>
<i>Issued on</i>	<b>2010-10-06</b>
<i>Validity date</i>	<b>2015-10-05</b>
<i>Approved by</i>	<b>Oscar M. Jordán (Managing Director)</b>

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