



Help with EOVE™-70 Settings



IPPB


The EOVE™-70 is a versatile device combining complementary treatments



Mode explanations

INEX

Mechanical in-exsufflation delivers positive pressure, followed very quickly by negative pressure, which assists coughing.




IPPB = Intermittent Positive Pressure Breathing


This *recruitment mode* delivers an air volume beyond the patient's current volume. It enables the patient to reach the inspiratory reserve volume. This improves pre-cough volume and mobilizes the thoracic wall.

The selected flowrate modulates the inspiratory cycle time

Two possible options depending on the treatment objective



The **high-flow** setting **effectively delivers** a volume of air over a reduced inspiratory time.



Slow insufflation allows **maximum intake** in all lung areas to reach the inspiratory reserve volume.

Indications and contraindications

Indications **INEX**

- Peak Cough Flow (PCF) when coughing < 160 L/min
 - Congestion
 - Thoracic wall diseases
 - Any patient with a PCF lower than 270 L/min will benefit from INEX
-

Indications **IPPB**

- Low PCR
- Congestion
- Decreased thoracic wall or pulmonary compliance
- Restrictive syndromes
- Forced vital capacity less than 60% of theoretical capacity

General contraindications applicable to both INEX and IPPB

Airway obstruction, hyperinflation, pneumothorax, severe emphysema and recent lobectomy, increased intracranial pressure, altered consciousness, lack of cooperation, severe bulbar palsy, cardiac instability, esophageal fistula, rib fractures.

Settings

INEX

Operating modes	Automatic/Manual	-	
Inspiratory pressure	5 to 70 cmH ₂ O	-	
Inspiratory time	0.5 to 5 sec	-	
Inspiratory slope	0 to 5	Speed at which the pressure will be reached: 1/5 th Ti, 2/5 th Ti, etc.	
Auto - Trigger	OFF/1 to 3	Allows the patient to trigger their cycle	
Pause	OFF/0.5 to 5 sec	Pause between expiration and inspiration	
PEEP	OFF/1 to 20 cmH ₂ O	Positive pressure during pause	
Expiratory pressure	0 to -70 cmH ₂ O	-	
Expiratory time	0.5 to 5 sec	-	
Oscillations	Inspiration/Expiration	Amplitudes	OFF/3 levels
		Frequency	OFF/4 to 20 Hz

Settings

IPPB

Inspiratory flowrate	5 to 100 L/min	Insufflation speed
Inspiratory flowrate	Continuous flowrate/ Decelerating flowrate	<i>Continuous flowrate</i> : constant flowrate during inspiration <i>Decelerating flowrate</i> : flowrate gradually decreases to the target pressure
Inspiratory trigger	OFF/8 levels	Allows the patient to trigger their cycle
Maximum pressure	10 to 50 cmH ₂ O	Pressure at which inspiration stops The patient can then exhale
Maximum time	0.5 to 20 sec	If the patient has not reached the max pressure at the end of the maximum inspiratory time, the cycle will proceed to expiration Maximum inspiratory phase time
Expiratory slope	5 levels	Gradual decrease in expiratory pressure to the set PEEP
PEEP	OFF/1 to 20 cmH ₂ O	Positive pressure during pause

To rinse the circuit, we recommend setting the expiratory slope to 1 and PEEP to 4 cmH₂O.

Starting a session: define comfort settings

IPPB

Support the patient to inhale a larger volume than their current volume.

Reassure the patient:

1

As you start breathing in through the mouthpiece, the EOVE™-70 will deliver air to you.

2

Allow yourself to be inflated by the device, **without forcing yourself to breathe in.**

3

When the device stops delivering air, you can exhale.

4


The exercise will take place at your pace: you decide when to start breathing in again.

duration: 5 to 15 min

Starting a session: setting comfort settings

INEX

Help the patient identify the inspiratory and expiratory phases so that they can synchronize as much as possible.



▶ Installations and recommendations for treatment

▶ Settle the patient into position based on their functional state:

- In a half-seated position in a bed
- Sitting in a chair

▶ Opt for a position where the practitioner can easily apply manual coughing assistance in addition to the instrumental assistance

If necessary, choose a personalized signal with the patient to indicate when to pause the insufflations (e.g. a wink)


▶ Ideally, use a mask that covers the nose and mouth with the INEX function

Guide the patient with clear instructions:

- Inhale or inflate
- Cough hard

► Objective

To simulate an effective cough to clear the patient's lungs



Settings

Recommended settings:

Installation settings	
Patient on NIV	Patient not on NIV
Initiate at: $P_i = PIP + 5 \text{ cmH}_2\text{O}$ $P_e = PIP + 10 \text{ cmH}_2\text{O}$	Initiate at: $P_i = 15 \text{ cmH}_2\text{O}$ $P_e = - 20 \text{ cmH}_2\text{O}$
$T_i = T_e = T_i (\text{NIV}) + 0.5 \text{ s}$ Pause = 1 sec	$T_i = T_e = 1 \text{ sec}$ Pause = 1 sec

Adjustment:

Inspiratory and expiratory pressures and times must be individualized and gradually increased until efficacy is achieved.

The T_i/T_e ratio must be adjusted according to the disease and the context.

The P_{insp} must be less than the P_{exp} .

Dose

Depends on the extent of the congestion,
example: 1 to 4 times a day; 2 to 5 cycles,
 pause when the patient needs a break.

Post-treatment clinical assessment

Improvement of the PCF, audible secretions,
 secretions in the upper airways.


Proximal clearance

Cough assist


► Objective

To increase inspiratory volume to facilitate effective coughing.


IPPB enhances pre-cough volume to facilitate coughing.



High-flow inspiration



Active deep inspiration



Secretion clearance

Settings

Recommended settings:

Flowrate: 30–50 L/min Max pressure: 30 to 40 cmH₂O

Adjustment:

Carried out by the physiotherapist during the session if there is an increase in bronchial hyperactivity:

- The flowrate and pressure must be individualized and gradually increased until a satisfactory pre-cough volume is achieved that enables effective coughing
- PEEP: in the absence of abdominal deficiencies, increase PEEP gradually

Dose

1 to 4 times a day depending on the congestion
Perform 2 to 5 movements followed by a cough;
pause when the patient needs a break

Post-treatment clinical assessment

Sputum, displayed volume increased, patient clinical improvement, auscultation.

Target patients


Congested patients with neuromuscular impairment: myopathy, tetraplegia.
Any congested patient whose dyspnea and/or exhaustion makes manual drainage ineffective.

Distal clearance

Mobilize secretions

► Objective

To mobilize secretions present in the distal parts of the lung areas.



Low-flow inspiration

Assisted active inspiration

Expiration: secretion movement

Settings

Recommended settings:

Flowrate: 20 to 30 L/min Max pressure: 25 to 35 cmH₂O

Adjustment:

If the pressure is reached too quickly or if the current volume displayed is too low: increase the inspiratory time by decreasing the flowrate and/or by adjusting the pressure gradually

Dose

1 to 4 times a day depending on the extent of the congestion.

Perform 2 to 5 insufflations followed by a cough; pause when the patient needs a break

Post-treatment clinical assessment

Sputum, displayed volume increased, patient clinical improvement, auscultation.

Target patients

Congested patients with neuromuscular impairment: myopathy, tetraplegia

Any congested patient whose dyspnea and/or exhaustion makes manual drainage ineffective.

Respiratory function optimization

Prevention of ventilation disorders

IPPB

► Objective

To maintain or develop ventilation in hypoventilated lung areas.



Slow inspiration

Assisted active slow inspiration

Posture to maximize inspiration

Settings

Recommended settings:

Flowrate: 10 to 30 L/min Max pressure: 25 to 35 cmH₂O

Adjustment:

Performed according to the mobility observed during auscultation of the hypoventilated area:

- Flowrate: minimum
- Maximum pressure: gradual titration of the pressure in order to optimize the inspired volume

Dose

10 to 30 minutes per day, 5 days a week

Target patients

- Neuromuscular diseases
- Multiple disabilities
- Pre-surgery, post-surgery
- Any patient whose respiratory function is severely decreased, with sustained volume reduction


Respiratory function optimization

Fight restrictions


IPPB

► Objective


To mobilize the chest cavity and let in as much air as possible.



Slow inspiration



Assisted active inspiration



Posture to maximize inspiration

Settings

Recommended settings:

Flowrate: 20 to 35 L/min Max pressure: 30 to 40 cmH₂O

Adjustment:

Depending on the volumes measured:

- Flowrate: increase or decrease depending on how the patient feels
- Max pressure: if resistance is high, do not hesitate to significantly decrease the pressure

Dose

10 to 30 minutes per day, 5 days a week

Target patients

- Neuromuscular diseases
- Multiple disabilities
- Pre-surgery, post-surgery
- Any patient whose respiratory function is severely decreased, with sustained volume reduction

References

- 1 Respiratoire, 2014. Comparison of three cough-augmentation techniques in neuromuscular patients: mechanical insufflation combined with manually assisted cough, insufflation-exsufflation alone and insufflation-exsufflation combined with manually assisted cough. Lacombe M1, Del Amo Castrillo L, Boré A, Chapeau D, Horvat E, Vaugier I, Lejaille M, Orlikowski D, Prigent H, Lafaso F.
- 2 AARC Clinical Practice Guideline. Intermittent Positive Pressure Breathing: 2003 Revision & Update. Sorenson HM, Shelledy DC, AARC.
- 3 Pediatric Pulmonology, 2006. IPPB-Assisted Coughing in Neuromuscular Disorders. Dohna-Schwake C, Ragette R, Teschler H, Voit T, Melies, U.
- 4 Arch. Bronconeumol., 2014. Comparison of intermittent positive pressure breathing and temporary positive expiratory pressure in patients with severe chronic obstructive pulmonary disease. Nicolini A1, Mollar E2, Grecchi B3, Landucci N2.

Contact

Air Liquide Medical Systems
Parc de Haute Technologie
6, rue Georges Besse
92182 Antony Cedex
France

www.device.airliquidehealthcare.com



Air Liquide Healthcare is a global leader in medical gases, home healthcare, hygiene products and specialty healthcare ingredients. Its mission is to provide its customers, throughout the hospital-to-home care pathway, with medical products, specialty ingredients and services that help protect vulnerable lives.

EOVE™-70 provides treatment for patients who are unable to clear their secretions by themselves. It offers an insufflation-exsufflation (INEX) mode and a pressure relieving mode (IPPB – Intermittent Positive Pressure Breathing) for adults and children – Class IIb medical device – CE 0459 – Manufacturer: EOVE™

Read the user manual carefully.

Air Liquide Medical Systems, Antony. Nanterre Trade and Companies Register No. 348 921 735