LANAFORM

HEALTH

WBPM-110

WRIST BLOOD
PRESSURE
MONITOR

2-60
MEMORIES

WHO
CLASSIFICATION

IRREGULAR HEART
RHYTHM ALERT

MEDICAL
DEVICE

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ENGLISH

Presentation

Thank you for choosing the WBPM-110 from Lanaform.

▲ Please read all the instructions before using your product, particularly these basic safety instructions.

☒ The photographs and other images of the product in this manual and on the packaging are as accurate as possible, but do not claim to be a perfect reproduction of the product.

Contents

- 01 Important information
- 02 Contents and display indicators
- 03 Intended use
- 04 Contraindication
- 05 Package contents
- 06 Product description.
- 07 Specifications
- 08 Notice
- 09 Setup and operating procedures
- 10 Maintenance
- 11 Explanation of symbols on unit
- 12 Limited warranty
- 13 Service center
- 14 Electromagnetic compatibility information
- 15 Advice regarding the disposal of waste

01. Important information

Normal Blood Pressure Fluctuation


All physical activity, excitement, stress, eating, drinking, smoking, body posture and many other activities or factors (including taking a blood pressure measurement) will influence blood pressure value. Because of this, it is mostly unusual to obtain identical multiple blood pressure readings.

Blood pressure fluctuates continually, day and night. The highest value usually appears in the daytime and lowest one usually at midnight. Typically, the value begins to increase at around 3:00AM, and reaches to highest level in the daytime while most people are awake and active.

Considering the above information, it is recommended that you measure your blood pressure at approximately the same time each day.

Too frequent measurements may cause injury due to blood flow interference, please always relax a minimum moment of 1 to 1.5 minutes between measurements to allow the blood circulation in your arm to recover. It is rare that you obtain identical blood pressure readings each time.

02. Contents and display indicators

- 01 Day and time display
- 02 Systolic
- 03 Diastolic
- 04 Pulse rate
- 05 Irregular heartbeat symbol
- 06 Low battery indicator
- 07 The pressure in cuff is instable, or much remnant air in cuff
- 08 Blood pressure level classification indicator
- 09 Blood pressure level classification indicator
- 10 **MEM** button
- 11 **START** button
- 12 LCD
- 13 Cuff

03. Intended Use

Fully Automatic Electronic Sphygmomanometer is for use by medical professionals or at home and is a non-invasive blood pressure measurement system intended to measure the diastolic and systolic blood pressures and pulse rate of an adult individual by using a non-invasive technique in which an inflatable cuff is wrapped around the wrist.

04. Contraindication

▲ It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer.

05. Package contents

- 01 Manual
- 02 Wrist blood pressure monitor
- 03 Battery 2×AAA

06. Product Description

Based on Oscillometric methodology and silicon integrated pressure sensor, blood pressure and pulse rate can be measured automatically and non-invasively. The LCD display will show blood pressure and pulse rate. The most recent 2×60 measurements can

be stored in the memory with date and time stamp. The Electronic Sphygmomanometer corresponds to the below standards: IEC 60601-1 Edition 3.1 2012-08/ EN 60601-1:2006/A1:2013 (Medical electrical equipment – Part 1: General requirements for basic safety and essential performance), IEC 60601-1-2:2014/ EN 60601-1-2:2015 (Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests), IEC 80601-2-30:2009+AMD1:2013/EN 80601 2-30:2010/A1:2015 (Medical electrical equipment – Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers) EN 1060-1: 1995 + A2: 2009 (Non-invasive sphygmomanometers - Part 1: General requirements), EN 1060-3: 1997 + A2: 2009 (Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro mechanical blood pressure measuring systems); ISO 81060-2:2013 (Non-Invasive Sphygmomanometers - Part 2: Clinical Validation Of Automated Measurement Type).

07. Specifications

- 01 Product name: WBPM-110
- 02 Model: LA090207
- 03 Classification: Internally powered, Type BF applied part, IP22, No AP or APG, Continuous operation
- 04 Machine size: Approx. 83mm × 74mm × 26mm (3 9/32" × 2 29/32" × 1 1/32")
- 05 Cuff circumference: 14cm-19.5cm (5 1/2"-7 11/16")
- 06 Weight: Approx. 71g (2 1/2oz.) (exclude batteries and cuff)
- 07 Measuring method: Oscillometric method, automatic inflation and measurement
- 08 Memory volume: 2×60 times with time and date stamp

- 09 Power source: batteries: 2 × 1.5V --- SIZE AAA
- 10 Measurement range:
 - a Cuff pressure: 0-300mmHg
 - b Systolic: 60-260mmHg
 - c Diastolic: 40-199mmHg
 - d Pulse rate: 40-180 beats/minute
- 11 Accuracy:
 - a Pressure: ±3mmHg
 - b Pulse rate: ±5%
- 12 Environmental temperature for operation: 10°C-40°C (50°F-104°F)
- 13 Environmental humidity for operation: ≤85%RH
- 14 Environmental temperature for storage and transport: -20°C-50°C (-4°F-122°F)
- 15 Environmental humidity for storage and transport: ≤85%RH
- 16 Environmental pressure: 80kPa-105kPa
- 17 Battery life: Approx. 200 times.
- 18 All components belonging to the pressure measuring system, including accessories: Pump, Valve, LCD, Cuff, Sensor

● These specifications are subject to change without notice.

08. ▲ Notice

- 01 Read all of the information in the operation guide and any other literature in the box before operating the unit.
- 02 Stay still, calm and rest for 5 minutes before blood pressure measurement.
- 03 The cuff should be placed at the same level as your heart.
- 04 During measurement, neither speak nor move your body and arm.
- 05 Measuring on left wrist for each measurement.
- 06 Please always relax at least 1 or 1.5 minutes between measurements to allow the blood circulation in your arm to recover. Prolonged over-inflation (cuff pressure exceed 300 mmHg

or maintained above 15 mmHg for longer than 3 minutes) of the bladder may cause ecchymoma of your arm.

- 07 Consult your physician if you have any doubt about below cases:
 - a The application of the cuff over a wound or inflammation diseases;
 - b The application of the cuff on any limb where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present;
 - c The application of the cuff on the arm on the side of a mastectomy;
 - d Simultaneously used with other monitoring medical equipments on the same limb;
 - e Need to check the blood circulation of the user.
- 08 ▲ This Electronic Sphygmomanometer is designed for adults and should never be used on infants or young children. Consult your physician or other health care professionals before use on older children.
- 09 Do not use this unit in a moving vehicle; this may result in erroneous measurement.
- 10 Blood pressure measurements determined by this monitor are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard Institute, Electronic or automated sphygmomanometers.
- 11 Information regarding potential electromagnetic or other interference between the blood pressure monitor and other devices together with advice regarding avoidance of such interference please see part ELECTROMAGNETIC COMPATIBILITY INFORMATION. It is suggested that the blood pressure monitor be kept 30 cm away from other wireless devices, such as WLAN unit, cell phone, microwave oven, etc. It can't be used near active HF SURGICAL EQUIPMENT and

the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high.

- 12** If Irregular Heartbeat (IHB) brought by common arrhythmias is detected in the procedure of blood pressure measurement, a  signal will be displayed. Under this condition, the Electronic Sphygmomanometers can keep function, but the results may not be accurate, it's suggested that you consult with your physician for accurate assessment. There are 2 conditions under which the signal of IHB will be displayed:
- a** The coefficient of variation (CV) of pulse period $>25\%$.
 - b** The difference of adjacent pulse period $\geq 0.14s$, and the number of such pulse takes more than 53 percentage of the total number of pulse.
- 13** Please do not use the cuff other than supplied by the manufacturer, otherwise it may bring biocompatible hazard and might result in measurement error.
- 14**  The monitor might not meet its performance specifications or cause safety hazard if stored or used outside the specified temperature and humidity ranges in specifications.
- 15**  Please do not share the cuff with other infected person to avoid cross-infection.
- 16** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
- a** Reorient or relocate the receiving antenna.
 - b** Increase the separation between the equipment and receiver.
 - c** Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
 - d** Consult the dealer or an experienced radio/TV technician for help.
- 17** Measurements are not possible in patients with a high frequency of arrhythmias.
- 18** The device is not intended for use on neonates, children or pregnant women. (Clinical testing has not been conducted on neonates, children or pregnant women.)
- 19** Motion, trembling, shivering may affect the measurement reading.
- 20** The device would not apply to the patients with poor peripheral circulation, noticeably low blood pressure, or low body temperature (there will be low blood flow to the measurement position).
- 21** The device would not apply to the patients who use an artificial heart and lung (there will be no pulse)
- 22** Consult your physician before using the device for any of the following conditions: common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, arterial sclerosis, poor perfusion, diabetes, pre-eclampsia, renal diseases.
- 23** The patient can be an intended operator.
- 24** This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

- 25 Attention that changes or modification not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.
- 26 Swallowing batteries and/or battery fluid can be extremely dangerous. Keep the batteries and the unit out of the reach of children and disabled persons.
- 27 If you are allergic to plastic/rubber, please don't use this device.

09. Setup And Operating Procedures

01 Battery Loading

- a Open battery cover at the back of the monitor.
- b Load two AAA size batteries. Please pay attention to polarity.
- c Close the battery cover.

When LCD shows battery symbol , replace all batteries with new ones.

Rechargeable batteries are not suitable for this monitor.

Remove the batteries if the monitor will not be used for a month or more to avoid relevant damage of battery leakage.


▲ Avoid the battery fluid to get in your eyes. If it should get in your eyes, immediately rinse with plenty of clean water and contact a physician.

 The monitor, the batteries and the cuff, must be disposed of according to local regulations at the end of their usage.

02 Clock And Date Adjustment

- a Once you install the battery or turn off the monitor, it will enter Clock Mode, and LCD will display time and date. See picture 2.

- b While the monitor is in Clock Mode, press the **START** and **MEM** button simultaneously for two seconds, then release the buttons, the monitor will be in Clock and Date Adjustment Mode.
- c In Clock and Date Adjustment Mode, the time format will blink at first. See picture 2-1. While time format is blinking, press the button **MEM** to change the time format. The default time format is 24h.
- d Press the button **START** repeatedly, the year (range is 2018-2099), month, day, hour and minute will blink in turn. See picture 2-2&2-3 & 2-4 & 2-5 & 2-6. While the number is blinking, press the button **MEM** to increase the number, keep on pressing the button **MEM**, the number will increase faster.
- e During adjusting clock and date, the monitor will go back to Clock Mode automatically when no button will be pressed within 60 seconds.
- f You can turn off the monitor by pressing **START** button when the minute is blinking, then the time and date is confirmed.

2-5

2-6

Note


- 01** The clock format could be set by user.
- 02** Table 2-1 instructs the conversion relations between 24 hour format and 12 hour format.

Table 2-1

24 HOUR FORMAT	12 HOUR FORMAT	24 HOUR FORMAT	12 HOUR FORMAT
0:00	12:00 AM	12:00	12:00 PM
1:00	1:00 AM	13:00	1:00PM
2:00	2:00 AM	14:00	2:00 PM
3:00	3:00 AM	15:00	3:00 PM
4:00	4:00 AM	16:00	4:00 PM
5:00	5:00 AM	17:00	5:00 PM
6:00	6:00 AM	18:00	6:00 PM
7:00	7:00 AM	19:00	7:00 PM
8:00	8:00 AM	20:00	8:00 PM
9:00	9:00 AM	21:00	9:00 PM
10:00	10:00 AM	22:00	10:00 PM
11:00	11:00 AM	23:00	11:00 PM

03 Connecting The Cuff To The Monitor

The cuff is attached to the monitor when it is packaged. Should the cuff become unattached, align the two plugs and four brackets of the cuff with the plug sockets and bracket sockets of the monitor and press the cuff to the monitor until the plugs and brackets are securely attached.

04 Applying The Cuff

- a** Place the cuff around a bare wrist 1-2cm above the wrist joint on the palm side of the wrist.
- b** While seated, place the arm with the cuffed wrist in front of your body on a desk or table with the palm up. If the cuff is correctly placed, you can read the LCD display.
- c** The cuff must be neither too tight nor too loose.

Note

- d** Please refer to the cuff circumference range in "SPECIFICATIONS" to make sure that the appropriate cuff is used.
- e** Measuring on left arm each time.
- f** Do not move your arm, body, or the monitor and do not move the rubber tube during measurement.
- g** Stay quiet, calm for 5 minutes before blood pressure measurement.
- h** Please keep the cuff clean. Clean the cuff by wet soft cloth and mild detergent if the cuff becomes dirty. Do not remove the cuff from the monitor. Clean the cuff after the usage of every 200 times is recommended.

Do not place the cuff around your arm if the arm has any inflammation, acute diseases, infections skin wounds.

05 Body Posture During Measurement

Sitting Comfortably Measurement

- Be seated with your feet flat on the floor, and don't cross your legs.
- Place palm upside in front of you on a flat surface such as a desk or table, with your elbow resting on a chair or table
- The middle of the cuff should be at the level of the right atrium of the heart.

06 Taking Your Blood Pressure Reading

- After applying the cuff and your body is in a comfortable position, press the **START** button. All display characters are shown for self-test. You can check the LCD display according to the picture. Please contact the service center if segment is missing. See picture 6-1.
- Then the current memory bank is blink. See picture 6-2. Press **MEM** button to change over to other bank. Confirm your selection by pressing **START** button. The current bank can also be confirmed automatically after 5 seconds with no operation.

- If the monitor has stored results, the LCD will momentarily display the most recent one. See picture 6-3. If no result has been stored, zero will appear on LCD. See picture 6-4.
- Then the monitor starts to seek zero pressure. See picture 6-5.
- Then the monitor inflates the cuff until sufficient pressure has built up for a measurement. Then the monitor slowly releases air from the cuff and carries out the measurement. Finally the blood pressure and pulse rate will be calculated and displayed on the LCD screen. The blood pressure classification indicator and Irregular heartbeat symbol (if any) will blink on the screen. See picture 6-6. The result will be automatically stored in the monitor.


- After measurement, the monitor will turn off automatically after 1 minute of no operation. Alternatively, you can press the **START** button to turn off the monitor manually.
- During measurement, you can press the **START** button to turn off the monitor manually.

ⓘ Please consult a health care professional for interpretation of pressure measurements.

07 Displaying Stored Results

- In Clock Mode, press **MEM** button, the monitor will display sign of current group. The amount of results in current user memory zone will be displayed. See picture 7-1. Press **START** button to switch group, press **MEM** to confirm current group. The current bank can also be confirmed automatically after 5 seconds with no operation.
- Then LCD will display the average value of all results in the current bank. See picture 7-2. If no result stored in the current user memory zone, LCD will display "0" for blood pressure and pulse rate. See picture 7-3.
- Press **MEM** button, the LCD will display the average value of the last three results in the current user memory zone. See picture 7-4. If no result stored in the current user memory zone, LCD will display "0" for blood pressure and pulse rate. See picture 7-5.
- Press **MEM** button again, the memory number in the current user memory zone will be displayed. See picture 7-6. If the monitor has no result stored in the current user memory zone, the LCD will display "0". See picture 7-7.
- Then the most recent result will be displayed with date and time stamp. See picture 7-8. Irregular heartbeat symbol (if any) and blood pressure classification indicator will blink at the same time. Press **MEM** button repeatedly to review the results measured previously. If the monitor has no result stored in the current user memory zone, the LCD will display "0" for blood pressure and pulse rate. See picture 7-9.
- When displaying the stored results, the monitor will turn off automatically after 1 minutes of no operation. You can also press the button **START** to turn off the monitor manually.

08 Deleting Measurements From The Memory

When any result is displaying (exclude average value displaying), keeping on pressing button **MEM** for three seconds, all results in the current bank will be deleted. See picture 8. Press the button **START** the monitor will be turned off.

09 Assessing High Blood Pressure For Adults

The following guidelines for assessing high blood pressure (without regard to age or gender) have been established by the World Health Organization (WHO). Please note that other factors (e.g. diabetes, obesity, smoking, etc.) need to be taken into consideration. Consult with your physician for accurate assessment, and never change your treatment by yourself.

Classification of blood pressure for adults

RANGE OF BLOOD PRESSURE VALUES	SYSTOLIC (IN MMHG)	DIASTOLIC (IN MMHG)	VALUE DISPLAYED
Level 3: acute hypertension	≥180	≥110	3
Level 2: intermediate hypertension	160-179	100-109	2
Level 1: mild hypertension	140-159	90-99	1
Normal high	130-139	85-89	N
Normal	120-129	80-84	N
Optimum	<120	<80	N

10 Technical Alarm Description

The monitor will show 'Hl' or 'Lo' as technical alarm on LCD with no delay if the determined blood pressure (systolic or diastolic) is outside the rated range specified in part SPECIFICACIONES. In this case, you should consult a physician or check if your operation violated the instructions.

The technical alarm condition (outside the rated range) is preset in the factory and cannot be adjusted or inactivated. This alarm condition is assigned as low priority according to IEC 60601-1-8.

The technical alarm is non-latching and need no reset. The signal displayed on LCD will disappear automatically after about 8 seconds.

Troubleshooting (1)

PROBLEM	POSSIBLE CAUSE	SOLUTION
LCD Display shows abnormal result	The cuff position was not correct or it was not properly tightened	Apply the cuff correctly and try again.
	Body posture was not correct during testing	Review the "BODY POSTURE DURING MEASUREMENT" sections of the instructions and re-test
	Speaking, arm or body movement, angry, excited or nervous during testing	Re-test when calm and without speaking or moving during the test
	Irregular heart-beat (arrhythmia)	It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer

Troubleshooting (2-1)

PROBLEM	POSSIBLE CAUSE	SOLUTION
LCD shows low battery symbol 	Low Battery	Change the batteries
LCD shows "ER 0"	Pressure system is unstable before measurement	Don't move and try again
LCD shows "ER 1"	Fail to detect systolic pressure	
LCD shows "ER 2"	Fail to detect diastolic pressure	

PROBLEM	POSSIBLE CAUSE	SOLUTION
LCD shows "ER 3"	Pneumatic system blocked during inflation	Connect the cuff correctly and try again. If the monitor is still abnormal, please contact the local distributor or the factory
LCD shows "ER 4"	Pneumatic system leakage during inflation	

Troubleshooting (2-2)


PROBLEM	POSSIBLE CAUSE	SOLUTION
LCD shows "ER 5"	Cuff pressure above 300mmHg	Measure again after five minutes. If the monitor is still abnormal, please contact the local distributor or the factory
LCD shows "ER 6"	More than 160 seconds with cuff pressure above 15 mmHg	
LCD shows "ER 7"	Inner Memory error	Take out batteries for five minutes, and then reinstall all batteries
LCD shows "ER 8"	Device parameter checking error	
LCD shows "ER A"	Pressure sensor parameter error	
No response when you press button or load battery	Incorrect operation or strong electromagnetic interference	

10. Maintenance

- 01 ▲** Do not drop this monitor or subject it to strong impact.
- 02 ▲** Avoid high temperature and solarization. Do not immerse the monitor in water as this will result in damage to the monitor.
- 03** If this monitor is stored near freezing, allow it to acclimate to room temperature before use.
- 04 ▲** The monitor requires 6 hours to warm from the minimum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is 20°C.
- 05** The monitor requires 6 hours to cool from the maximum storage temperature between uses until the the monitor is ready for its INTENDED USE when the ambient temperature is 20°C.
- 06 ▲** Do not attempt to disassemble this monitor.
- 07** If you do not use the monitor for a long time, please remove the batteries.
- 08** It is recommended the performance should be checked every 2 years or after repair. Please contact the service center.
- 09** Clean the monitor with a dry, soft cloth or a soft cloth squeezed well after moistened with water, diluted disinfectant alcohol, or diluted detergent.
- 10** No component can be maintained by user in the monitor. The circuit diagrams, component part lists, descriptions, calibration instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated repairably can be supplied.
- 11** The monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or three years, and the cuff integrity is maintained after 1,000 open-close cycles of the closure.
- 12** It is recommended the cuff should be disinfected 2 times every week if needed (For example, in hospital or in clinique). Wipe the inner side (the side contacts skin) of the cuff by a soft cloth squeezed after moistened with Ethyl alcohol (75-90%), then dry the cuff by airing.
- 13** Not servicing/maintenance while the monitor is in use.

11. Explanation of Symbols On Unit

The following symbols appear in these instructions for use:

SYMBOL	DESCRIPTION
	Warning instruction indicating a risk of injury or damage to health.
	Safety note indicating possible damage to the unit/accessory.
	Note on important information.
	Symbol for "THE OPERATION GUIDE MUST BE READ" (The sign background colour: blue. The sign graphical symbol: white).
	Symbol for "SERIAL NUMBER"
	Symbol for "TYPE BF APPLIED PARTS" (The cuff is type BF applied part)
	Symbol for "ENVIRONMENT PROTECTION - Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice".
	Symbol for "COMPLIES WITH MDD93/42/EEC REQUIREMENTS"
	Manufacturer
	Symbol for "EUROPEAN REPRESENTATION"
	Level of protection against water penetration.
	Date of manufacture.

12. Limited warranty

LANAFORM guarantees this product against any material or manufacturing defect for a period of two years from the date of purchase, except in the circumstances described below.

The LANAFORM guarantee does not cover damage caused as a result of normal wear to this product. In addition, the guarantee covering this LANAFORM product does not cover damage caused by abusive or inappropriate or incorrect use, accidents, the use of unauthorized accessories, changes made to the product or any other circumstance, of whatever sort, that is outside LANAFORM's control.

LANAFORM may not be held liable for any type of circumstantial, indirect or specific damage.

All implicit guarantees relating to the suitability of the product are limited to a period of two years from the initial date of purchase as long as a copy of the proof of purchase can be supplied.

On receipt, LANAFORM will repair or replace your appliance at its discretion and will return it to you. The guarantee is only valid through LANAFORM's Service Centre. Any attempt to maintain this product by a person other than LANAFORM's Service Centre will render this guarantee void.

13. Service Center

 Andon Health Co., Ltd.
No. 3 Jinping Street, YaAn Road,
Nankai District, Tianjin 300190, China
Tel: +86 22 87611660

 iHealthLabs Europe SAS
36 Rue de Ponthieu, 75008, Paris, France
Tel: +33 1 42661559

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14. Electromagnetic Compatibility Information

Table 1 – Emission

PHENOMENON	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT
RF emissions	CISPR 11 Group 1 Class B	Home Healthcare environment
Harmonic distortion	IEC 61000-3-2 Compliance	Home Healthcare environment
Voltage fluctuations and flicker	IEC 61000-3-3 Compliance	Home Healthcare environment

Table 2 – Enclosure Port

PHENOMENON	BASIC EMC STANDARD	IMMUNITY TEST LEVELS
		Home Healthcare environment
Electrostatic Discharge	IEC 61000-4-2	± 8kV contact ± 2kV, ±4kV, ±8kV, ±15 kV air
Radiated RF EM filed	IEC 61000-4-3	10V/m 80MHz-2.7GHz 80% AM at 1kHz
Proximity filed from RF wireless communications equipment	IEC 61000-4-3	Refer to table 3
Rated power frequency magnetic fields	IEC 61000-4-8	30A/m 50Hz or 60Hz

Table 3 – Proximity filed from RF wireless communications equipment

TEST FREQUENCY (MHZ)	BAND (MHZ)	IMMUNITY TEST LEVELS
		Professional Healthcare facility environment
385	380-390	Pulse modulation 18Hz, 27V/m
450	430-470	FM, ±5kHz deviation, 1kHz sine, 28V/m
710	704-787	Pulse modulation 217Hz, 9V/m
745		
780		
810	800-960	Pulse modulation 18Hz, 28V/m
870		
930		
1720	1700-1990	Pulse modulation 217Hz, 28V/m
1845		
1970		
2450	2400-2570	Pulse modulation 217Hz, 28V/m
5240	5100-5800	Pulse modulation 217Hz, 9V/m
5500		
5785		

15 Advice regarding the disposal of waste

All the packaging is composed of materials that pose no hazard for the environment and which can be disposed of at your local sorting centre to be used as secondary raw materials. The cardboard may be disposed of in a paper recycling bin. The packaging film must be taken to your local sorting and recycling centre.

When you have finished using the device, please dispose of it in an environmentally friendly way and in accordance with the law.

Prior to disposal, please remove the battery and dispose of it at a collection point so that it can be recycled.

Under no circumstances should spent batteries be disposed of with household waste.

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Ltd., No. 3 Jinping
Europe SAS, 36 Rue