

Pipep Kit



- Polycarbonate valve, ABS connectors and resistors
- Latex free, DEHP and BPA free
- The kit is composed by one two way valve: IN= inspiration OD 22mm OUT= expiration OD 22mm ID 15mm, and a mouthpiece which is connectable to a mask or to a tracheostomy tube. Contained in the kit are 8 resistors for PEP or PIP with color code with the following diameters:
 1,5 mm black,
 2,0 mm white,
 2,5 mm yellow,
 3,0 mm blue,
 3,5 mm green,
 4,0 mm orange,
 4,5 mm pink,
 5,0 mm brown;
- Manometer connector, OD 15 mm
- Mask not included in the kit
- Single patient ISO

Code	Description	Items per box
KM-803	Pipep kit	24

Crystal Masks to use with Pipep Kit



Code	Description	Size	Color Code
KM-200	Mask	Ultra large - 7	Trasparent
KM-201	Mask	Large - 6	Red
KM-202	Mask	Medium - 5	Yellow
KM-203	Mask	Small - 4	Green
KM-204	Mask	Pediatric- 3	White
KM-205	Mask	Infant - 2	
KM-206	Mask	Neonatal - 1	
KM-207	Mask	Pre-neonatal -0	



Manometer

Manometer for monitoring the operating pressures

Code	Description	Measure	Peso	Items per box
57301537	Manometer	- 50 + 50 cmH ₂ O	100 g	1
57303499	Manometer	- 30 + 30 cmH ₂ O	100 g	1
57355227	Manometer	- 100 + 150 cmH ₂ O	100 g	1

Tri-Gym

Lung exercises



- Polycarbonate-made, PVC balls, DEHP and BFA free, polyethylene tube and mouthpiece
- Latex free
- For inspiration and expiration
- Two adjustable valves controlling independently inspiration and expiration
- Three measuring cylinders for the flow control
- Flow control is done through the three colored balls, which indicate the amount of flow per second 600 cc green, blue 900 cc, 1200 cc red.
- Corrugated tube and mouthpiece
- Single patient use

Code	Description	Item per box
KM-805	Tri-Gym	24