



## **AirZone** Peak Flow Meter Technical specifications

**Tradition and Innovation** - Visionary thinking and a fascination with technology have guided us throughout our history to develop innovative products of outstanding reliability. Anticipating trends to improve the quality of life.

**HS** CLEMENT CLARKE  
INTERNATIONAL

# AirZone PFM

## Technical specifications



### AirZone PFM

Dimensions	139 x 44 x 35mm	
Weight	44g	
Measurement Range	60 - 800 L/min (EU Scale)	60 - 720 L/min (ATS Scale)
Accuracy	±10% or 10 L/min	
Repeatability	<5 L/min	

Product	Part Codes
AirZone PFM (EU Scale) in a box	3130160
AirZone PFM (EU Scale) in a poly bag	3131001
AirZone PFM (ATS Scale) in a box	3130119
AirZone PFM (ATS Scale) in a poly bag	3130057
AirZone Windmill Trainer (box of 50)	3104262

## Product information

The AirZone PFM is a pocket size peak flow meter with the same accuracy as a Mini-Wright Standard, with the added value of the zoning function for easy asthma management. Suitable for both adults and children. There is a Windmill Trainer available ideal for encouraging young children.

- Available with ATS and EU Scales
- Individually calibrated for absolute accuracy
- Conforms to NAEP guidelines
- Zoning system allows self management
- Windmill Trainer available, a helpful incentive for young children