

Wrist blood pressure monitor BW 310

Art.-No. 51070 | EAN 40 15588 51070 0

Accurate blood pressure measurement on the wrist

60 memory spaces for each of 2 users

Display items: systolic, diastolic, pulse, date, time



Application | Features

Accurate blood pressure measurement on the wrist

Display items: systolic, diastolic, pulse, date, time

Measurements are classified according to the WHO* evaluation system using a traffic-light colour scale Arrhythmia indicator

60 memory spaces for each of 2 users

Calculates average value of data of the last 3 results

Certified medical device**

Technical Data

Size: approx. 7,2 x 7,5 x 3 cm

Weight: 150 g

Measuring method: Oscillometric

Measuring range:

- Blood pressure: 30 – 250 mmHg
- Pulse: 40 – 180 beats/min.

Batteries: 2 x 1.5 V (AAA)

Wrist cuff for 13,5–21,5 cm wrist circumferences

Package Content

Wrist blood pressure monitor

Batteries

Storage case

Logistical Data

Sales packaging

1 piece | tbd | tbd

Packaging unit

tbd | tbd | tbd

* World Health Organisation

** Certification according to MDD (Medical Device Directive)