

Upper arm blood pressure monitor MTC

Art.-No. 51139 | EAN 40 15588 51139 4

Accurate blood pressure measurement on the upper arm

60 memory slots

Voice output in selectable languages: DE, GB, ES, NL, FR, IT, PT

Extra large, easy-to-read display



Application | Features

Accurate blood pressure measurement on the upper arm

Display items: systolic, diastolic, pulse, date, time

Extra large, easy-to-read display

Measurements are classified according to the WHO® evaluation system using a traffic-light colour scale

60 memory slots

Voice output in selectable languages: DE, GB, ES, NL, FR, IT, PT

Automatic power off after 1 minute

Adjustable sound volume

Optional mains or battery operation (mains adapter optional available)

With storage pouch

Certified medical device**

* World Health Organisation

**Certification according to MDD (Medical Device Directive)

Technical Data

Size: approx. 10.8 x 15.1 x 6.3 cm

Weight: approx. 0.30 kg

Measuring method: Oscillometric

Measurement range:

Blood pressure: 0 – 300 mmHg

Pulse: 30 – 180 beats/min.

Batteries: 4 x 1.5 V (AA)

Scope of Delivery

Upper arm blood pressure monitor

Cuff for 22 – 32 cm upper arm circumferences

Batteries

Storage bag

Logistical Data

Sales packaging

1 piece | 32 x 22 x 8 cm | 1.2 kg

Packaging unit

4 pieces | 36.5 x 25 x 37 cm | tbd