

Upper arm blood pressure monitor MTS

Art.-No. 51152 | EAN 40 15588 51152 3

Accurate blood pressure measurement on the upper arm

60 memory slots for each of 2 users

Arrhythmia indicator

Extra large, easy-to-read display



Application | Features

Accurate blood pressure measurement on the upper arm

Display items: systolic, diastolic, pulse, date, time

Extra large, easy-to-read display

Arrhythmia indicator

Measurements are classified according to the WHO* evaluation system using a traffic-light colour scale

60 memory slots for each of 2 users

Automatic power off after 1 minute

Optional mains or battery operation (mains adapter optional available)

With storage pouch

Certified medical device**

* World Health Organisation

** Certification according to MDD (Medical Device Directive)

Technical Data

Size: approx. 19.5 x 22.5 x 26.5 cm

Weight: approx. 0.30 kg

Measuring method: Oscillometric

Measurement range:

Blood pressure: 0 – 300 mmHg

Pulse: 30 – 180 beats/min.

Batteries: 4 x 1.5 V (AA)

Scope of Delivery

Upper arm blood pressure monitor

Cuff for 22 – 30 cm upper arm circumferences

Batteries

Storage bag

Logistical Data

Sales packaging

1 piece | 17.7 x 11.8 x 8.8 cm | 0.3 kg

Packaging unit

4 pieces | 22.7 x 20.1 x 27.5 cm | 3.6 kg