

Upper arm blood pressure monitor BU 510

Art.-No. 51160 | EAN 40 15588 51160 8
PZN 9888613

Accurate blood pressure measurement on the upper arm

Measurements are classified according to the WHO*

Very large cuff

Very large digits for easy reading



Application | Features

Accurate blood pressure measurement on the upper arm

Display items: systolic, diastolic, pulse, date, time

Arrhythmia indicator

Measurements are classified according to the WHO* evaluation system using a traffic-light colour scale

90 memory slots for each of 2 users

Calculates average value of all stored data

Very large cuff (arm circumferences 22 – 36 cm)

Very large digits for easy reading (up to 22 mm in height)

With storage bag

Certified medical device**

* World Health Organisation

**Certification according to MDD (Medical Device Directive)

Technical Data

Size: approx. 13 x 11 x 6 cm

Weight: approx. 0.30 kg

Measuring method: Oscillometric

Measurement range:

Blood pressure: 0 – 299 mmHg

Pulse: 40 – 199 beats/min.

Batteries: 4 x 1.5 V (AA)

Scope of Delivery

Upper arm blood pressure monitor

Cuff for 22 – 36 cm upper arm circumferences

Batteries

Storage bag

Logistical Data

Sales packaging

1 piece | 16 x 8.5 x 13 cm | 0.6 kg

Packaging unit

4 pieces | 34.8 x 16.4 x 13.4 cm | 2.5 kg