

## Wrist blood pressure monitor HGV

Art.-No. 51220 | EAN 40 15588 51220 9

Accurate blood pressure measurement on the upper arm

120 memory slots for each of 2 users

Voice output in selectable languages: DE, GB, ES, NL, FR, IT, PT

Extra large, easy-to-read display



### Application | Features

Accurate blood pressure measurement on the wrist

Display items: systolic, diastolic, pulse, date, time

Extra large, easy-to-read display

Measurements are classified according to the WHO\* evaluation system using a traffic-light colour scale

Arrhythmia indicator

120 memory slots for each of 2 users

Calculates average value of all stored data

Voice output in selectable languages: DE, GB, ES, NL, FR, IT, PT

Adjustable sound volume

With storage case

Certified medical device\*\*

\* World Health Organisation

\*\* Certification according to MDD (Medical Device Directive)

### Technical Data

Size: approx. 8.5 x 8.2 x 3.6 cm

Weight: approx. 0.13 kg

Measuring method: Oscillometric

Measurement range:

Blood pressure: 0 – 300 mmHg

Pulse: 0 – 180 beats/min.

Batteries: 2 x 1.5 V (AAA)

Circumference of cuff for wrist: 14–19.5 cm

### Scope of Delivery

Wrist blood pressure monitor · Batteries

Storage case

### Logistical Data

#### Sales packaging

1 piece | 8.5 x 8.2 x 3.6 cm | 0.75 kg

#### Packaging unit

4 pieces | 40 x 17.5 x 24.5 cm | tbd