

# MEDISANA®

## Inhalator IN 550 Pro

Art.-No. 54530 | EAN 40 15588 54530 6

For targeted treatment of diseases of the upper and lower respiratory tract

Highly efficient due to very precise nebulisation

Inhalation with mask, mouthpiece or nose adapter – ideal for children

Extra long hose – 2 m



extra long hose – 2m



comprehensive set of accessories



integrated accessory case



integrated cable storage



Compact and portable

### Application | Features

Nebulisation with compressed air technology

For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)

Highly efficient due to very precise nebulisation

Particularly effective thanks to high proportion of respirable aerosols

Integrated accessory case

Comprehensive set of accessories:  
Inhalation with mask, mouthpiece or nose adapter – ideal for children

Extra long hose – 2 m

Certified medical device\*

\*Certification according to MDD (Medical Device Directive)

### Technical Data

Size: approx. 33 x 20 x 13.8 cm

Weight: approx. 2 kg

Power output: 60 watts

Power supply: 230V/50 Hz

Filling volume inhalant: max. 10 ml

Nebulizing output: ≤ 0,2 ml/min

Particle size: >60% ≤ 5 µm

Hose length: 2 m

### Package Content

Inhalator

Mouth piece

Masks for children and adults

Connecting piece

Nose adapter for children and adults

5 replacement filters

2 replacement fuses

### Logistical Data

Sales packaging

1 piece | tbd | tbd

Packaging unit

4 pieces | tbd | tbd