

The MetaNeb™ 4.0 System

The MetaNeb™ System is designed for secretion mobilisation, lung expansion and aerosol delivery in a single therapy for the treatment and prevention of pulmonary atelectasis.

Technical Specifications

Controller

weight	2.2 kg
height (mounted to stand)	25.4 cm
width (mounted to stand)	38.1 cm
diameter (mounted to stand)	10.2 cm

Bracket Arm*	0.5 kg
---------------------	--------

Stand

weight	6 kg
rolling diameter	63.5 cm
base and pole	127.64 cm

Oxygen Supply Hose

colour	green
length	25.4 cm
fittings on each end	DISS

Single Patient Use (SPU) Circuit**


colour	blue
length	152.4 cm

Power Requirements

oxygen source	3.45 bar
---------------	----------

* Attaches to standard 2.5 cm diameter pole. Height adjustable throughout the pole length.

** Disposable Circuit materials are not manufactured with natural rubber latex. Product available in single and 5-pack.



Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

MetaNeb™ is a registered trademark of Comedica, Inc. This document is destined solely for use by healthcare professionals. Medical devices shown in this brochure are intended for use with patients in departments of healthcare establishments.

These medical devices are regulated health products which, pursuant to such regulation, bear a CE mark. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN165306-01, 08 February 2016

Enhancing outcomes for patients and their caregivers:

Hill-Rom