

Infrared-Multifunctional thermometer TM 750

Art.-No. 76140 | EAN 40 15588 76140 9

Precise measurement of the body temperature in the ear and on the forehead in seconds

Additional function: measurement of ambient temperature, liquids and surfaces

Visual fever alarm

Automatic storage of the last 30 measurement results



Visual fever alarm
– green LED: temperature normal , red LED: raised temperature



Application | Features

Precise measurement of the body temperature in the ear and on the forehead in seconds

Additional function: measurement of ambient temperature, liquids and surfaces (0 – 100° C)

Visual fever alarm – green: temperature normal, red: raised temperature

Acoustic signal at the beginning and end of temperature reading

Automatic storage of the last 30 measurement results

Date and time

Certified medical device*

Technical Data

Size: approx. 13.4 x 4.2 x 4.7 cm

Weight: approx. 65kg

Measurement range object: 0-100 °C

Measurement range body temperature: 34-43 °C

Scaling: 0.1°C

Accuracy of measurement:

- Forehead: +/- 0.3°C from 34 - 43 °C

- Ear: +/- 0.2°C from 35.5 - 42 °C

- Object: ≤30 +/-1.5 °C
>30 +/-5%

Batteries: 2 x 1.5 V (AAA)

Scope of Delivery

Infrared-Multifunctional thermometer

Batteries

Storage bag

Logistical Data

Sales packaging

1 piece | 10.8 x 16.6 x 6 cm | tbd

Packaging unit

tbd | tbd | tbd

*Certification according to MDD (Medical Device Directive)