

## Thermometer TM 700


Art.-No. 77040 | EAN 40 15588 77040 1

Precise measurement of the body temperature:  
oral, axillary, rectal

Flexible tip

Measurement result in 10 seconds

Watertight



### Application | Features

Precise measurement of the body temperature:  
oral, axillary, rectal

Acoustic signal at the beginning and end of  
temperature reading

Acoustic fever alarm

Automatic storage of the last memory result

Watertight

Automatic power off after 10 minutes

Digital LCD display

Flexible tip

Conversion from °C to °F

Measurement result in 10 seconds

Easy battery exchange

With storage case

Certified medical device\*

\*Certification according to MDD (Medical Device Directive)

### Technical Data

Size: approx. 13.9 x 2.2 x 1.3 cm

Weight: approx. 0.12 kg

Measurement range: 32,0 – 42,9°C

Scale increment: 0,1°C

Battery: 1 x 1,5 V (AG3/LR41)

### Scope of Delivery

Thermometer

1 micro battery (AG3/LR41)

Storage case

### Logistical Data

#### Sales packaging

1 piece | tbd | tbd

#### Packaging unit

tbd | tbd | tbd