

PURE NEB®

CLASSIC



High performance nebulizing system AERODJINN®

Treatment of broncho-pulmonary pathologies
(BPCO, DDB, asthma, bronchitis, pneumopathy, cystic fibrosis)

This compact device, combined with ExtraNEB+®, allows nebulizing all medicines types (pneumatic compressor and continuous aerosol production – MMAD: 3.9 µm, with ExtraNEB+® and AERODJINN®+ compressor).

Performances with ExtraNEB+®:

2 ml = 4 min. - 5 ml = 11 min.,

MMAD: 3.9 µm,

Inhalable mass: 25%.

Characteristics:

Weight: 1.9 kg


Dimensions: 22x12x17.5 cm

Warranty: 3 years

BRONCHO-PULMONARY
TREATMENT



Combined with ExtraNEB+® nebulizing kit, AERODJINN® compressor allows a uniform deposition in lungs.



Scintigraphy on a healthy patient

Latex and phthalate free accessories:



EXTRA NEB+®, with mouth piece



Paediatric maks (RP) included

Optional accessories:



Infant nebulizing kit NL9MBB



Adult mask, R



Medical Device class IIa. Read carefully the instructions for use delivered with the medical device.
Documentation for Health Professionals. Non contractual pictures.

A-COMPS.EN.009.0715

AEROSOLTHERAPY BY DTF MEDICAL



Diffusion Technique Française
19, rue de la Presse - CS 60132
42003 Saint-Etienne Cedex 1 - France
Tel. : (33) 4 77 74 51 11 - Fax : (33) 4 77 79 67 72
E-mail : dtf@dtf.fr - www.dtf.fr

