

AEROcheck[®]

Fingertip Pulse Oximeter



Pulse Oximetry

Spot-Check measurement | SpO₂ and pulse

AEROcheck®

Fingertip Pulse Oximeter

AEROcheck® - P011-FCH Fingertip Pulse Oximeter

Fingertip Pulse Oximeter for measurement of partial oxygen saturation.

Can be used anywhere: comfortable, quick and precise.

Exact

- Display:
SpO2 | Pulse frequency | Pulse intensity
Plethysmogram | Battery capacity
- Clinically approved measurement accuracy

Clear

- Dual colored display with 360 degree adjustments by pressing button (6 display modes)

Handy

- Portable and lightweight (58 g)
- Comfortable transport with attached storage bag

Persistent

- Approx. 30 hours of continuous operation
- Optical low-battery alarm

Durable

- High quality material and finishing
- Finger case made of high quality silicone
- Stable battery case with hinged cover



Versatile field of application

- Hospital / Clinic
- Homecare
- Ambulance
- Other medical facilities
- Private
- Sports

Used for: Pediatric and adult

„Unique in its' quality!“



Soft-Slip finger case



Robust spring mechanism



Hinged battery cover

Simple usage

Automatic start/stop function

3 Steps - 1 Result

- Open
- Put finger in
- Check measurements

Done!



One-Button operation
for display settings
and brightness



360°
adjustable
display



Comfortable haptics
through silicone finger case

Storage and Transport

Flexible | Versatile | Protected

Storage case

- Carrying string
- Belt loop
- Fixation eyelet

„Always and
everywhere
with you!“



Tether



Belt loop



High quality and durable nylon material

Scope of delivery and Technical Data

- ① Fingertip Pulse Oximeter | ② Batteries | ③ Carrying string
- ④ Storage case | ⑤ Instructions for use

Specification	Model	
	P011-FCH (REF: HP011-FCH)	
Dimensions (WxHxD)	60 x 34 x 38 mm	
Weight	60 g incl. batteries	
Internal power supply	Batteries: 2 pieces AAA - alkaline 1.5 V	
Operating time	approx. 30 h	
SpO ₂	Display	0 % - 99 %
	Calibration	70 % - 99 %
	Resolution	1 %
Puls	Display	0 bpm - 254 bpm
	Calibration	30 bpm - 235 bpm
	Resolution	1 bpm
Display	SpO ₂ , Plethysmogram, Pulse frequency, Pulse display	
Alarm (optical)	Low-battery display	
Operation temperature	+5° C - +40° C	
Storage temperature	-20° C - +55° C	
Relative humidity	5 % - 93 %	
Atmospheric pressure	86 kPa - 106 kPa	



Scope of Delivery

Fingertip Pulse Oximeter

AEROcheck®

P011-FCH

Item-No.	Description	Dimensions (WxHxD)	Weight
HP011-FCH	AEROcheck® - P011-FCH, Fingertip Pulse Oximeter for children and adults, dual color OLED-Display with 60 x 34 x 38 mm 58 g (incl. batteries) 6 Display mode and automatic switch off function, incl. Nylon safety case as well as 2 x AAA alkaline batteries	60 x 34 x 38 mm	60 g (incl. batteries)

HUM

HUM Gesellschaft für Homecare
und Medizintechnik mbH
Zum Pier 79 | D-44536 Lünen
FON +49 (0)2 31/ 88 08 85-0 | FAX +49 (0)2 31/ 88 08 85-58
Mail: sales@hum-online.de | http://www.hum-online.de

AEROcheck® - Fingertip pulse oximeter
PP-P011-FCH_en | Status 23.11.2015 | Revision 2
Changes and errors reserved. Illustrations are approximate.
© HUM Gesellschaft für Homecare und Medizintechnik mbH