



Swedish Hospital Supply AB

BA-TUBE Product Datasheet

Issuer:	Date:	Edition:	Basic UDI-DI:	Designation:
SHS, Swedish Hospital Supply, AB	2022-04-07.	2.	B-17350026250201	BA-TUBE

Product description

BA-TUBE (Bronchitis Aid Tube) is a - PEP (Positive Expiratory Pressure) and PIP (Positive Inspiratory Pressure) device – an all-in-one function consisting of a tube with an adjustable top for setting the correct resistance.

BA-TUBE has been developed in collaboration with Pulmonary Physicians and Physiotherapists at University Hospital US in Linköping and KI Karolinska University Hospital in Stockholm.

Intended use

BA-TUBE is intended to train and strengthen the respiratory system of the user. BA-TUBE is used by patients suffering from COPD (Chronic Obstructive Pulmonary Disease) and is used as a pre- and postoperative respiratory aid.

Characteristics and benefits

BA-TUBE is a flexible and discreet breathing trainer that facilitates patient compliance with the caregiver's recommended exercise program. BA-TUBE's discreet design makes it easy to use anywhere in society and when the need arises without attracting attention. With such benefits, BA-TUBE increases the patient's chances of having a good therapeutic effect and compliance with the caregiver's training program.

BA-TUBE requires no loose or extra parts to adjust the resistance, due to the mounted top with 8 different resistors.

BA-TUBE has a long service life and in the long term a low daily cost.

BA-TUBE, with its two components (Tube and Top) is easy to use, easy to disinfect and meets the hygiene requirements that hospitals place on a breathing aid.

The TUBE

Material: Polypropylene, White, Green text printing

Manufacturing: Injection molded



Swedish Hospital Supply AB

TOP adjuster

Material: Polypropylene, Green
Manufacturing: Injection molded

Stick/Plunger

Material: Polypropylene, black
Manufacturing: Injection molded

Packaging

Article nr	UDI-DI- Cardboard Box	Number/carton	Weight/carton
100003-A	17350026250027	1 st BA-TUBE	15 g
	17350026250034	50 pcs/box	approx. 1,5 kg
	17350026250041	500 pcs/box	approx. 15 kg

Material Cardboard: Unbleached corrugated board.

Label on carton: UDI-DI for each packaging volume, LOT, date of manufacture (year, month, day), expiration date (year, month, day).

Production

The product is produced in accordance with **ISO 9001**. Quality tests are carried out routinely.

Country of origin

Sweden

Sterile product

No

Shelf life

10 years, unbroken packaging.

1 year warranty from the date of purchase regarding material defects.



Swedish Hospital Supply AB

Storage

Store dry at room temperature.

Handling Waste

Recycling hard plastic product

Other information

<https://swehospitalsupply.se/en/>



Secure-Innovative-Sympathetic