

Body composition analyzer BS 475

Art.-No. 40432 | EAN 40 15588 40432 0

Measurement of weight, body fat, body water, muscle and bone mass

Convenient „step-on“ switch

Tread area of 6 mm safety glass

Maximal load capacity: 180 kg



Flat design – only 22 mm



Application | Features

Measurement of weight, body fat, body water, muscle and bone mass

Determination of the BMI value

Convenient „step-on“ switch

Tread area of 6 mm safety glass

Large easy-to-read LCD display with white luminous digits

Toggles between kg, lb and st

Four high-precision strain gauge sensors for accurate results

10 memory slots

Automatic switch-off

Overload indicator

Maximal load capacity: 180 kg

Technical Data

Size: approx. 30 x 30 x 2.2 cm

Weight: approx. 1.6 kg

Measuring range up to 180 kg, 396 lb, 28 st 4 lb

Graduation: 100 g, 0,2 lb, 1/4 lb

Display size: 83 x 45 mm

Batteries: 2x 3V (CR2032)

Scope of Delivery

Body composition analyzer

Batteries

Logistical Data

Sales packaging

1 piece | 32.8 x 32.5 x 4.4 cm | tbd

Packaging unit

4 pieces | 34,3 x 19,1 x 34,5cm | 7,2 kg