

smiths

# Graseby® MS Syringe Drivers

when accuracy matters



INFUSION MANAGEMENT

# GRASEBY® SYRINGE DRIVERS



MS16A Syringe Driver rate set in mm/hr

<b>Typical Infusion Range</b>	30 mins – 60 hrs
<b>Variable Speed</b>	0 – 99 mm/hr



MS26 Syringe Driver rate set in mm/day


<b>Typical Infusion Range</b>	12 hours – 60 days
<b>Variable Speed</b>	0 – 99 mm/day
<b>Bolus</b>	0.23 mm of travel

Famed for their simplicity and reliability, the Graseby® MS Drivers have been the choice of healthcare professionals in hospitals, nursing homes, palliative care and community settings for over 25 years.

## Applications

Suitable for both IV and subcutaneous infusion, typical drug infusion include: analgesics, sedatives, diuretics, bronchodilators, chemotherapy agents.

- a. Overlaps to prevent objects levering case apart.
- b. Clear window to view syringe pusher.
- c. Clear window to view syringe barrel.
- d. Syringe Luer protected.
- e. Slot for Infusion Set.



- f. Key-operated lock for security.
- g. Locating 'pip' to secure pump.
- h. Hard polycarbonate (shatterproof) material.
- i. Separating joint to aid cleaning.
- j. Clear area for viewing set rate and operating instructions (on rear of pump).



## Syringe Driver Packed Set

Syringe driver, shoulder holster, syringe cover, rate adjuster, battery and instruction manual.

Product	Product Code
MS16A	0105-0504
MS26	0113-0001

Note: These part-numbers refer to the UK variant.

## The Graseby® MS Driver Lock Box

- Consists of hard polycarbonate (shatterproof) material incorporating clear viewing windows for monitoring the syringe and rate settings.
- Lightweight and portable. Stable on a flat surface or fits into the Graseby® MS Syringe Driver shoulder holster if required.
- Key lock access offers complete control and increased confidence to the authorised clinician.

## Accessories

Product	Product Code
100cm Infusion Set PVC line with 25 gauge needle	0105-0029
Subcutaneous Infusion Pack BD 10ml syringe, Tegaderm dressing, alcoholic wipe transfer needle and 100cm infusion set	0105-0117
Syringe driver – non slip base	0105-0108
Cover	0105-0529
Holster Washable, soft fabric holster with cover	0105-0027
Instruction Manual	0105-0549
Rate Adjuster Tool to turn slotted rate switches	0113-0023
Training Pack Full package of presentation, instruction and testing material	TPF-00130
Me and My Syringe Driver (Patient Information Booklet)	0060-1012
Poster Memory aid for Syringe Driver set up	TPF-00302
Lock Box Additional protection from fluid ingress. Fits the MS Driver Shoulder Holster and allows stable operation on a flat surface	0105-0640

## Fact:

The Archer Fish, which live in the Mangroves of South East Asia, are famed for their ability to shoot a highly accurate jet of water from their mouths in order to knock down prey.

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.


For further information please call your local Smiths Medical distributor or Smiths Medical on +44 (0)1303 260551

# smiths

**Smiths Medical International Ltd**

Hythe, Kent CT21 6JL UK  
Tel: +44 (0)1303 260551 Fax: +44 (0)1303 266761  
www.smiths-medical.com

Smiths, Graseby are registered trademarks of Smiths Group plc.  
Literature No. LIT/MD2536



# smiths

Smiths Medical - a part of Smiths Group plc