


# ***Alpha 300***

The reference IPPB



**NEW**



[www.airliquidemedicalsystems.com](http://www.airliquidemedicalsystems.com)



# Alpha 300

The reference IPPB  
Pediatric and Adult

Alpha 300 is a multi-purpose device to be used by physiotherapists for the prevention of respiratory disorders. It delivers positive pressure during inspiration and switches over to expiration as soon as the pre-set pressure has been reached.

## 1 Optimization of respiratory function

Surpass the patient's maximum inspiration

➡ Increase his vital capacity

## 2 Help in pulmonary recruitment

Introduce large volumes of air to recruit the badly or non-ventilated parts of the lung

➡ Improve haematosis and pulmonary expansion

## 3 Assistance with lung clearing

Deliver to the patient a volume of air large enough to:


- Mobilize the deepest secretions
- Cough or sigh to clear out these secretions

➡ Allow an effective patient expiration (increase of sputum and reduction of drainage effort)

## Alpha 300 highlights

- User friendly
- Simple and complete follow-up on LCD screen
- Can be used in hospitals, in the physiotherapists offices and at home

➡ Alpha 300 is a unique tool for an autonomous respiratory rehabilitation





## Target pathologies

- Neuromuscular diseases and other restrictive illnesses
- COPD, mucoviscidosis
- Pre and post surgery therapies

### References:

KC039600 Alpha 300

Breathing circuit:

KG020900	Circuit adulte avec nébuliseur 1 m Ø 22 mm, x20	
KG021000	Circuit adulte sans nébuliseur 1 m Ø 22 mm, x20	



**AIR LIQUIDE**

[www.airliquide.com](http://www.airliquide.com)