



**PARI PEP® System I/II**



**PARI PEP® S System**



# PARI PEP® System

## Time-saving combination of physiotherapy and inhalation treatment

**The key features at a glance**

- Mobilises secretions and thus facilitates clearing of the lungs
- Stabilises the respiratory tract by adjustable expiration resistance: counteracts bronchial collapse
- Pressure display for setting expiration resistance and monitoring during PEP treatment (PARI PEP® System II)
- Applicable in physiotherapy without nebuliser

<b>PARI PEP® S System</b>	<b>PARI PEP® System I</b>	<b>PARI PEP® System II</b>
The PARI PEP® S System can be used with all PARI nebulisers that are currently available commercially (except the PARI LC SPRINT BABY, the PARI BABY and the PARI TIA)	For <b>all</b> PARI Nebulisers (except PARI LC SPRINT® family, PARI BABY and PARI TIA) Without manometer	For <b>all</b> PARI Nebulisers (except PARI LC SPRINT® family, PARI BABY and PARI TIA) With manometer
<b>Item No. 018G4000</b>	<b>Item No. 018G6101</b>	<b>Item No. 018G6201</b>

**Accessories for the PARI PEP® System**



**Manometer**  
0–100 mbar, including Pressure Hose  
**Item No. 018B1600**

**PARI PEP® Tubing adapter**  
Pack of 3, made of silicon, with one adapter for manometer  
**Item No. 041G4580**

