



NEW

G3/G5 Cart

Product Specifications

Catalog Number and Price

Instructions for Use

oxygen. anytime. anywhere.®

INO GEN ONE® 3



INO GEN ONE® 5



PRODUCT SPECIFICATIONS



oxygen. anytime. anywhere.®




Specifications

WEIGHT


3.25 lbs (1.47kg)



17.68in
(44.90cm)
COLLAPSED



34in
(86.36cm)
PARTIAL HEIGHT



41in
(104.14cm)
EXTENDED HEIGHT

5.45in
(13.84cm)

12.51in
(31.77cm)

oxygen. anytime. anywhere.®

CATALOG #

oxygen. anytime. anywhere.®



Catalog # and Prices

Catalog Number	Description
CA-501	Inogen G3/G5 Cart



oxygen. anytime. anywhere.®



INSTRUCTIONS FOR USE

oxygen. anytime. anywhere.®



Instructions for Use with G3 Device

INOGEN ONE® G3



Slide the G3 Carry Bag handle over the cart handle.



Intake, exhaust vents and cannula barb should be facing forward.



Large rubber-based wheels to allow smooth gliding up and down stairs and along rough uneven pavement and surfaces.

G3 should seat firmly in lower base of cart to avoid from falling out.

oxygen. anytime. anywhere.®

inogen

Instructions for Use with G5 Device

INOGEN ONE® G5



Opening on back of cart for exhaust ventilation.

Slide the G5 Carry Bag over the cart handle through the sleeve on back of carry bag.



Front opening flap of carry bag and cannula barb should be facing forward.



Intake vents

- > located on each side of the G5
- > visible through side mesh of carry bag.

G5 should seat firmly in lower base of cart to avoid from falling out.

oxygen. anytime. anywhere.®



