

Catalogue Number

Disposable range	
KM222100	Hospital NV Full Face Disposable mask - Size S
KM222200	Hospital NV Full Face Disposable mask - Size M
KM222300	Hospital NV Full Face Disposable mask - Size L
Reusable range	
KM218800	Hospital NV Full Face Reusable mask - Size S
KM218900	Hospital NV Full Face Reusable mask - Size M
KM219000	Hospital NV Full Face Reusable mask - Size L
KM222400	Disposable headgear X10

Compatibility with ventilators

	Monnal T50	Monnal T60	Monnal T75	Extend	Other valve ventilators
					
Compatibility	✓	✓	✓	✓	✓



Contact

Air Liquide Medical Systems - Parc de Haute Technologie
 6, rue Georges Besse - 92182 Antony cedex
 Tel.: +33 (0)1 40 96 66 00 - Fax: +33 (0)1 40 96 67 00



Hospital NV Full Face Disposable mask

The solution to benefit from our comprehensive expertise in non-invasive ventilation



Safety

Efficiency


Time-saving

- Provides non-invasive, comfortable ventilation
- Guarantees complete safety for the user
- Saves time by simplifying traceability

YN152111 - V1.1 - 29/01/2013 - Raffaello Fischer


www.airliquidemedicalsystems.com

Hospital NV Full Face Disposable mask



- Mask available in 3 sizes
- Sizing-gauge attached on the packaging
- Package includes mask, cushion, and headgear assemblies

Fitting recommendations



Duration of use

- The authorised duration of use is 7 days.
- There is no approved disinfection or sterilisation method for this medical device.
- For use on multiple patients, please use our range of reusable masks.

