

# In-Check DIAL

**Many patients do not use their inhalers correctly, failing to inhale at the correct flow rate. Guidelines call for inhaler review. In-Check Dial enables the healthcare professional to coach patients to use their inhalers properly.**

## Product Information

The In-Check Dial measures peak Inspiratory flow rate. It can be set to simulate the resistance characteristics of many common inhalers.

Patients can be trained, using In-Check Dial, to inhale at a flow rate known to be suitable for a particular type of inhaler. The In-Check Dial shows the flow rate achieved alongside a scale showing the clinically effective flow ranges for different types of inhalers.

## What is it Used for?

The In-Check Dial is used for inhaler review to coach patients to improve their inhaler technique.



## Suitability

Hospital ✓

Clinic ✓

## Product Features

- Individually calibrated to ensure a high degree of accuracy
- Easy to clean by simple washing
- Resistance selector simulates common inhalers
- Use with one way mouthpieces
- Describes correct technique for each inhaler and clinically effective range

## In-Check DIAL

### Technical Specifications

Measurement Range: 15 -120 L/min

Accuracy: 10% or  $\pm 10$  L/min

Repeatability: 5 L/min

### Part Codes

In-Check DIAL without Information Scale: **3109900**

In-Check DIAL with Information Scale: **3109910**

In-Check DIAL 5 T&S: **3109969**

In-Check Inhaler Assessment Kit: **3109760**


### Two-Way Mouthpieces

Cost effective mouthpieces that allow the In-Check DIAL to be used as a multi-patient device. Cleaning is recommended between each patient when using these mouthpieces.

### Part Codes

Adult Disposable Two Way Mouthpieces

Box of 100 (5 x 20): **3125030**

Box of 500: **3122003**

Paediatric Disposable Two Way Mouthpieces

Box of 40: **3125005**

Box of 500: **3123030**

In-Check Paediatric Adaptor for use with Paediatric Disposable Two-Way Mouthpieces (Pack of 1): **3125027**

## Consumables



### One-Way Inspiratory Mouthpieces

Allows In-Check DIAL to be used as a multi-patient device. One Way valve helps to prevent cross infection.

### Part Codes

Adult Disposable Inspiratory One-Way Valve Mouthpieces

Box of 100 (5 x 20): **3122061**

Adult Disposable Inspiratory One-Way Valve Mouthpieces

Box of 200: **3122063**

Adult Disposable Inspiratory One-Way Valve Mouthpieces

Box of 500: **3122065**

Paediatric Mouthpiece Adaptor (Pack of 1): **3122049**

Paediatric Mouthpiece Adaptor (Pack of 20): **3122120**