

Coach[®] 2 and CliniFLO[®] Incentive Spirometers Improving Quality of Life



RESPIRATORY CARE



Coach® 2 Volume Incentive Spirometers

Coach® 2 makes incentive spirometry easier for patients and staff alike. The one-way valve ensures patients inhale, rather than exhale into the unit. A highly-visible piston, easy to understand graphics (indicating correct inspiratory flowrate), and instructions in the base, help patients perform and monitor their own post-surgical breathing exercises without direct supervision. Other features include an O₂ connection for supplemental oxygen, convenient handle, flexible popple tubing and bedrail holder.

Product Benefits

- Ensures patients inhale, rather than exhale into the unit.
- Easily adjustable for each patient's use.
- Can be seen by patients emerging from the effects of anaesthesia.
- Easy to train.
- Stays with the patient for maximum compliance.
- Accommodates patients requiring supplemental oxygen.
- Saves space.
- Colourful deep-sea characters, games, puzzles and stickers.

Coach® 2 Incentive Spirometers

ORDERING INFORMATION

PART No.	PRODUCT DESCRIPTION	Volume
22-4000	Coach®2 One way valve	4000ml
22-4001	Coach®2	4000ml
22-2500	Coach®2 One way valve	2500ml
22-2501	Coach®2	2500ml
22-2000	Coach®2 Kids One way valve	2000ml



CliniFLO® Low-Flow Incentive Spirometers

CliniFLO® is ideal for therapy in geriatric, paediatric or weakened patients. With flow settings from 100ml/sec to 600ml/sec, virtually any patient can sustain the minimum inspiratory effort required for effective therapy. Slow inspirations enhance collateral ventilation and minimise patient discomfort when performing post-surgical breathing exercises. The O₂ port makes it easy to provide oxygen during therapy. ClinifLO® unit comes with preassembled popple tubing, mouthpiece and instructions.

Product Benefits

- Can accommodate virtually any patient.
- Reduces the possibility of contamination.
- Provides immediate patient feedback.
- Easy to train.
- Reduces the chance that the setting will be changed inadvertently.

CliniFLO® Low-Flow Incentive Spirometers

ORDERING INFORMATION

PART No.	PRODUCT DESCRIPTION
22-1200	ClinifLO®

For further information please visit:
www.smiths-medical.com/respiratorycare

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

Smiths Medical International Ltd

Hythe, Kent CT21 6JL UK
 Tel: +44 (0)1303 260551 Fax: +44 (0)1303 266761

Smiths Medical ASD, Inc.

160 Weymouth Street
 Rockland, MA 02370 (USA)
 Tel: (800) 848-1757 Fax: (800) 621-2950
www.smiths-medical.com

Smiths Medical, part of the global technology business Smiths Group

The Smiths Medical and Portex design marks; Coach and ClinifLO are trademarks of the Smiths Medical family of companies. The symbol © indicates the trademark is registered in the U.S. Patent and trademark office and certain other countries.
 Literature No. LIT/AM3005

smiths medical