

Blood glucose monitor MediTouch®

Art.-No. 79025 | EAN 40 15588 79025 6

PZN 5358261

Measuring time within 5 seconds

0.6 µL blood

360 memory slots

Simple menu navigation using symbols



**360
MEMORY**



**3 YEARS
WARRANTY**



Complete starter set in convenient pouch



Application | Features

Measuring time within 5 seconds

0.6 µL blood

Ultra slim design

Large LCD display (32 x 42 mm)

360 memory slots

Calculates average values for 7/14/30/60/90 days

Simple menu navigation using symbols

Pre and post meal markers

Certified medical device*

Technical Data

Size: approx. 7.5 x 4.4 x 1 cm

Weight: approx. 40 g

Battery: 1 x 3 V (CR 2032)

Scope of Delivery

Blood glucose monitor

10 test strips | Lancing device

10 lancets | Control solution

AST cap | 5 alcohol pads

Pouch | User manual

Quick reference guide | Diabetic diary

Battery

Logistical Data

Sales packaging

1 piece | 17.7 x 7.3 x 11.8 cm | 0.24 kg

Packaging unit

4 pieces | 30.5 x 18.7 x 13.4 cm | 1.6 kg

*Certification according to IVD Directive respectively MDD (Medical Device Directive)