

AEROTube® - BSS12-IPPB | Universal IPPB tubing system

IPPB tubing system for a respiratory therapy with intermittent positive pressure and simultaneous medicament nebulisation. Compatible with popular pressure inhalation devices*.

Efficient IPPB respiratory therapy in hospitals or at home!

22 mm PE corrugated tube

- Light weight
- Flexible
- Safe

Noise-optimized patient valve

- 360° rotatable expiration opening
- For a gentle therapy

Nebulizer (20 ml)

- Efficient aerosol production

2 regulation hoses

- Expiration
- Medicament nebulisation

Anatomy accessible mouthpiece



Newly developed patient valve „whisperingly quiet!“



Easy to use

- Quick Setup
- Comfortable handling
- Lightweight PE (160 g)

Economical and reliable

- Inexpensive
- No reprocessing costs

Transparent

- Monitoring of the complete system

Optimal respiratory therapy

- Constant volume respiration and simultaneous aerosol inhalation with less effort
- Production of dense aerosol
- Very high patient comfort



* not suitable for devices in which the adjustment of the exhalation resistance is made on stenosis caps