

Ombra* Table Top Compressor, 230 Volt – Euro and G Plug Models

*Designed for optimal performance with
AeroEclipse* XL BAN* Nebulizer*



Silicone cover keeps the power switch free of dust and dirt

Contoured handle with nebulizer cup holder for added convenience

Carrying case keeps treatment accessories organized



Manufacturer: Trudell Medical International

Origin: Canada

Ordering Information:

	Euro Plug	G Plug
Part Number	10551193006	10551593006
Single Unit GTIN	(01) 1 0762860 05511 4	(01) 3 0762860 05515 6
Carton GTIN	(01) 2 0762860 05511 1	(01) 4 0762860 05515 3
Quantity in Carton	6	6
Packaging	Retail carton	Retail carton
Package Contents	Ombra* 230 Volt Table Top Compressor with Euro Plug, two AeroEclipse* XL BAN* Nebulizers with mouthpiece and tubing, five replacement air filters, carrying case, and instructions	Ombra* 230 Volt Table Top Compressor with G Plug, two AeroEclipse* XL BAN* Nebulizers with mouthpiece and tubing, five replacement air filters, carrying case, and instructions
Product Reference Number (REF)	105511	105515

Device Characteristics:‡

- Electrical Requirements: 230 Volt AC
- Power: 180 VA Maximum
- Frequency: 50 Hz
- Weight: 1.59 kg (3.5 lb)
- Dimensions: 17.8 cm X 15.5 cm X 10 cm (7 in X 6.1 in X 3.9 in)
- Sound Level: 60 dB(A)
- Power Cord Length: 1.52 m (5 ft)
- Duty Cycle: 30 minutes on/30 minutes off
- Equipment Classification: Class II
- Operating Temperature Range: +5°C to +40°C (+41°F to +104°F)
- Operating Humidity Range: 10 to 90% (non-condensing)
- Operating Atmospheric Range: 75 kPa (11 psi) to 101.3 kPa (15 psi)
- Storage/Transport Temperature Range: -25°C to +70°C (-13°F to +158°F)
- Storage/Transport Humidity Range: 10 to 90% (non-condensing)
- Minimum Pressure: 103 kPa (15 psi)
- Nominal Operating Pressure with the **AeroEclipse* XL BAN*** Nebulizer: 140 kPa (20 psi)
- Nominal Operating Flow with the **AeroEclipse* XL BAN*** Nebulizer: 4L/min
- Free Flow Air Rate: 10 L/min
- Operating Life: 1,095 hours
- CE marked

Cleaning Method:‡

- Wipe exterior of housing with a damp cloth

Disinfection Methods:‡

- Wipe the case, switch, and handle with a tissue dipped in one of the following:
 - 70% isopropyl alcohol
 - Bleach solution
 - 3M[†] SoluPrep[†] Solution
- Wipe the case, switch, and handle with a PDI[†] Super Sani-Cloth[†] Germicidal Disposable Wipe

Maintenance:‡

- Change the air filter each year

‡ Refer to the [instructions for use](#) for additional information.



MD-474A-0924 * trade marks and registered trade marks of Trudell Medical International (TMI). † trade marks and registered trade marks of their respective companies. © TMI 2024. All rights reserved.

Manufactured by:



TRUDELL MEDICAL
INTERNATIONAL

+1-519-455-4862
info@trudellmed.com trudellmed.com

LEARN
MORE



345-GLO-EN