

AEROcheck®

Fingertip Pulse Oximeter



new
HUM - NOVELTY

Pulse Oximetry

Spot-Check measurement | SpO2 and pulse

AEROcheck® - Fingertip Pulse Oximeter

AEROcheck® - FP011-FCH: Handy Fingertip Pulse Oximeter for measurement of partial oxygen saturation.
Can be used anywhere: comfortable, quick and precise.

Exact

- Display: SpO₂ | Pulse frequency | Pulse intensity | Plethysmogram | Battery capacity
- Clinically approved measurement accuracy

Clear

- Dual coloured display with 360 degree adjustments by pressing button (6 display modes)

Handy

- Portable and lightweight (58 g)
- Comfortable transport with attached storage bag

Persistent

- approx. 30 hours of continuous operation
- Optical low-battery alarm

Durable

- High quality material and finishing
- Finger case made of high quality silicone
- Stable battery case with hinged cover



Used for:
Pediatric and adult



Versatile field of application

- Hospital / Clinic
- Homecare
- Other medical facilities
- Private
- Ambulance
- Sports

„Unique in its' quality!“



Soft-Slip finger case



Robust spring mechanism



Hinged battery cover

AEROcheck® - simple usage

3 steps - 1 result

Automatic On / Off function

- Open
- Put finger in
- Check measurements

Done!

One-Button operation for display settings and brightness



AEROcheck® - Storage and Transport

Flexible | Versatile | Protected

Storage case

- Carrying string
- Belt loop
- Fixation eyelet

„Always and everywhere with you!“



High quality and durable nylon material

AEROcheck® - Product Range and Technical Data

① Fingertip Pulse Oximeter | ② Batteries | ③ Carrying string | ④ Storage case | ⑤ Instructions for use



Specific.		Model	P011-FCH (REF: HP011-FCH)
Dimensions (WxHxD)			60 x 32 x 35 mm
Weight			58 g incl. batteries
Internal power supply			Batteries: 2 pieces AAA - alkaline 1,5 V
Operating time			approx. 30 h
SpO ₂	Display		0% - 99%
	Calibration		70% - 99%
Puls	Display		0 bpm - 254 bpm
	Calibration		30 bpm - 235 bpm
	Resolution		1 bpm
Display			SpO ₂ , Plethysmogram, Pulse frequency, Pulse display
Alarm (optical)			Low-battery display
Operation temperature			+5° C - +40° C
Storage temperature			-20° C - +55° C
Relative humidity			5 % - 93 %
Atmospheric pressure			86 kPa - 106 kPa

Product overview - Fingertip Pulse Oximeter

AEROcheck® - P011-FCH

Item-No.	Description	Dimensions (WxHxD)	Weight
HP011-FCH	Fingertip Pulse Oximeter for children and adults, dual colour OLED-Display with 6 Display mode and automatic switch off function, incl. Nylon safety case as well as 2 x AAA alkaline batteries	60 x 34 x 38 mm	58 g (incl. batteries)



AEROcheck®

HUM

HUM Gesellschaft für Homecare
und Medizintechnik mbH

Zum Pier 79 | D-44536 Lünen

FON +49 (0)2 31/ 88 08 85-0 | FAX +49 (0)2 31/ 88 08 85-58

Mail: sales@hum-online.de | <http://www.hum-online.de>

AEROCHECK® - Fingertip Pulse Oximeter
PP-P011-FCH_GB | Stand 26.09.2013 | Revision 1

© HUM Gesellschaft für Homecare und Medizintechnik mbH