

# AEROresc®

EASY Collar | Cervical Collar



# AEROresc®

## EASY Collar | Cervical Collar

### AEROresc® - EASY Collar | Cervical Collar

The AEROresc® - EASY Collar is a size adjustable cervical collar intended for safe immobilization of the cervical spine. It reliably restricts the freedom of movement of the cervical spine and takes over its stabilizing function.

The fine grid of both models assures easy and individual size adjustment even beyond standard sizes. Due to its self-raising chin support the cervical collar could be given a flat design and can easily be stored with minimum space requirements.

With only 2 models for any patient type (from infant up to tall adults) you are perfectly equipped for emergencies.

EASY Collar - Highest stability at maximum wearing comfort.

### Persuasive Advantages

- reliable restriction of freedom of movement
- lock-mechanism assures high stability
- symmetrical size adjustment due to guide rail
- easy and quick size selection due to finger technique
- self-raising anatomical chin support
- spacious access front and rearside for pulse control, airway management, view at larynx, palpation
- left and right-side fastening hook for mounting of O<sub>2</sub> tubing
- radioluscent, suitable for CT & MRI
- space saving storage due to flat design
- high quality material - hypallergenic, latex- and PHT-free
- instruction for use printed onto product - always at hand

### The Models



AEROresc® - EASY Collar Mini

Infants | Children | Small Adults

Dimensions: approx. 45 x 14 cm | Weight: approx. 70 g



AEROresc® - EASY Collar Adult

Adults: Neckless | Short | Regular | Tall

Dimensions: approx. 57 x 18 cm | Weight: approx. 120 g

2 versions: Mini | Adult



AEROresc® - EASY Collar Mini




AEROresc® - EASY Collar Adult

### Area of Application

- emergency service
- hospital

### Easy Handling



### Fastening Hooks for O<sub>2</sub> Tubing

- left- and right-side
- safe fixation of O<sub>2</sub> tubing

### Flexible Size Adjustment

- symmetrical size adjustment
- 16 level adult-version:  
1-4 neckless | 5-8 short | 9-12 regular  
13-16 tall
- 12 level mini-version:  
1-4 infants | 5-8 children  
9-12 small adults

### Comfortable Chin Support

- tilts up automatically during application
- anatomical and well-cushioned


### Lock System

- assures high stability
- striking buttons for easy detection

### Spacious access front- and backside





- palpation of cervical spine
- better air circulation



- pulse control
- larynx view
- airway management

## Scope of Delivery EASY Collar | Cervical Collar

### AEROresc® - EASY Collar | Cervical Collar

Product Code	Description		Dimensions / Weight
HIMT11-EC-K	AEROresc® - EASY Collar Mini, cervical collar for infants, children and small adults, SU = 20 pcs.		45 x 14 cm / 70 g
HIMT11-EC-E	AEROresc® - EASY Collar Adult, cervical collar for adults, SU = 20 pcs.		57 x 18 cm / 120 g

