

Respireo SOFT Baby

Supporting his first breath



Developed in collaboration with neonatologists and pediatric health experts, Respireo SOFT Baby is a baby dedicated mask up to approximately 15 months.


Babies faces differ greatly in their shapes from those of adults, in particular with a prominent forehead, no defined nose bridge and a growing & fragile bone structure. All the mask features are designed to address such specificities.

Entirely made of soft silicone, the cushion perfectly fits the facial characteristics of the patient and **ensures minimal pressure points**.

With 5 fixation points, the headgear is specifically designed to guarantee **masks stability and perfect seal** around the nose without tightening.

The **mask can be customized** without forehead support to adapt to the most demanding faces (malformation, highly convex foreheads etc.).

Respireo **SOFT** Baby



5 points headgear

Large straps and the addition of a 5th fixation strap on the forehead ensures homogenous bearing points and stability whatever the pressures or baby's position.



Silent and diffuse vent system

To fit into the peaceful environment of the baby.

Soft and comfortable cushion

With a unique and thin lip to guarantee both air tightness and minimal pressure on the face.

External port

To adapt to all kinds of non invasive treatments needs (measure tidal CO₂, oxygen supply etc.)

Rotary and flexible tubing at 360°

For a total sense of freedom.

References

Respireo SOFT Baby XS	KM225300
Respireo SOFT Baby S	KM225400
Respireo SOFT Baby headgear XS	KM225500
Respireo SOFT Baby headgear S	KM225600



Please read user manual.
Medical device class IIa, CE 0051.
Manufactured by Air Liquide Medical Systems Spa. Distributed by Air Liquide Medical Systems SA.