



Suction Catheters

A comprehensive family of easy to use atraumatic suction catheters for a range of procedures.



COVIDIEN

positive results for life™

Gentle to patients

A wide range of atraumatic suction catheters with choice of tip designs and optional connectors for different suction sources. Constructed with a frosted surface to simplify insertion and to an effective level of hardness to resist kinking.

Suction Catheter with Moulded Connector and Bevelled Tip

- Frosted surface for easy insertion
- Appropriate degree of hardness provides kink resistance
- Atraumatic, soft rounded, bevelled open tip
- Two lateral eyes with smooth edges
- Colour coded connector
- Effective length: 50cm

All products latex free

No.	Colour	Code
1.	06 ch Light Green	190 0601 1
2.	08 ch Blue	190 0801 1
3.	10 ch Black	190 1001 1
4.	12 ch White	190 1201 1
5.	14 ch Green	190 1401 1
6.	16 ch Orange	190 1601 1
7.	18 ch Red	190 1801 1
8.	20 ch Yellow	190 2001 1




Bevelled tip

Easier identification

Universal size colour coding for tips and connectors across the range.

BS EN ISO 8836; 2009



Suction Catheter Rounded Tip

All products latex free

- Frosted surface for easy insertion
- Appropriate degree of hardness provides kink resistance
- Ideal tip reducing the risk of tissue damage
- Atraumatic, soft rounded, open tip
- Two small relieving eyes
- Colour coded connector
- Effective length: 50cm




No.	Colour	Code
9.	● 06 ch Light Green	190 0685 1
10.	● 08 ch Blue	190 0885 1
11.	● 10 ch Black	190 1085 1
12.	● 12 ch White	190 1285 1
13.	● 14 ch Green	190 1485 1
14.	● 16 ch Orange	190 1685 1
15.	● 18 ch Red	190 1885 1
16.	● 20 ch Yellow	190 2085 1



Rounded tip

Suction Catheter with Vacuum Control and Bevelled Tip

All products latex free



- Frosted surface for easy insertion
- Appropriate degree of hardness provides kink resistance
- Atraumatic, soft rounded, bevelled open tip
- Two lateral eyes with smooth edges
- Colour coded vacuum control connector
- Effective length: 50cm

No.	Colour	Code
17.	● 06 ch Light Green	191 0611 1
18.	● 08 ch Blue	191 0811 1
19.	● 10 ch Black	191 1011 1
20.	○ 12 ch White	191 1211 1
21.	● 14 ch Green	191 1411 1
22.	● 16 ch Orange	191 1611 1
23.	● 18 ch Red	191 1811 1

Bevelled tip

Suction Catheter with Vacuum Control and Rounded Tip


- Frosted surface for easy insertion
- Appropriate degree of hardness provides kink resistance
- Metric type has accurate positioning depth marks at each cm starting 4cm from the distal end
- Ideal tip reducing the risk of tissue damage
- Atraumatic, soft rounded, open tip
- Two small relieving eyes
- Vacuum control connector
 - Transparent for easy observation of secretions
 - Colour coded
 - No contamination of finger
- Effective length: 50cm

All products latex free

No.	Colour	Code
24.	● 05 ch Grey	190 0551 1
25.	● 06 ch Light Green	190 0651 1
26.	● 07 ch Violet	190 0751 1
27.	● 08 ch Blue	190 0851 1
28.	● 10 ch Black	190 1051 1
29.	○ 12 ch White	190 1251 1
30.	● 14 ch Green	190 1451 1
31.	● 16 ch Orange	190 1651 1
32.	● 18 ch Red	190 1851 1



Rounded tip



Suction Catheter Coude Tip

All products latex free



- Angled tip for endobronchial suction procedures, providing easier access to the right or left bronchus
- Frosted surface for easy insertion
- Appropriate degree of hardness provides kink resistance
- Ideal tip reducing the risk of tissue damage
- Atraumatic, soft rounded, open tip
- Two small relieving eyes
- Marker on funnel connector indicating direction of angled tip
- Colour coded connector
- Effective length: 50cm

No.	Colour	Code
33.	● 10 ch Black	190 1065 1
34.	○ 12 ch White	190 1265 1
35.	● 14 ch Green	190 1465 1
36.	● 16 ch Orange	190 1665 1



Ordering Information

Suction Catheter with Moulded Connector and Bevelled Tip

No.	Colour	Code	Sales unit
1.	● 06 ch Light Green	190 0601 1	100
2.	● 08 ch Blue	190 0801 1	100
3.	● 10 ch Black	190 1001 1	100
4.	○ 12 ch White	190 1201 1	100
5.	● 14 ch Green	190 1401 1	100
6.	● 16 ch Orange	190 1601 1	100
7.	● 18 ch Red	190 1801 1	100
8.	● 20 ch Yellow	190 2001 1	100

Suction Catheter Rounded Tip

No.	Colour	Code	Sales unit
9.	● 06 ch Light Green	190 0685 1	100
10.	● 08 ch Blue	190 0885 1	100
11.	● 10 ch Black	190 1085 1	100
12.	○ 12 ch White	190 1285 1	100
13.	● 14 ch Green	190 1485 1	100
14.	● 16 ch Orange	190 1685 1	100
15.	● 18 ch Red	190 1885 1	100
16.	● 20 ch Yellow	190 2085 1	100

Suction Catheter with Vacuum Control and Bevelled Tip

No.	Colour	Code	Sales unit
17.	● 06 ch Light Green	191 0611 1	100
18.	● 08 ch Blue	191 0811 1	100
19.	● 10 ch Black	191 1011 1	100
20.	○ 12 ch White	191 1211 1	100
21.	● 14 ch Green	191 1411 1	100
22.	● 16 ch Orange	191 1611 1	100
23.	● 18 ch Red	191 1811 1	100






Suction Catheter with Vacuum Control and Rounded Tip

No.	Colour	Code	Sales unit
24.	● 05 ch Grey	190 0551 1	100
25.	● 06 ch Light Green	190 0651 1	100
26.	● 07 ch Violet	190 0751 1	100
27.	● 08 ch Blue	190 0851 1	100
28.	● 10 ch Black	190 1051 1	100
29.	○ 12 ch White	190 1251 1	100
30.	● 14 ch Green	190 1451 1	100
31.	● 16 ch Orange	190 1651 1	100
32.	● 18 ch Red	190 1851 1	100






Suction Catheter Coude Tip

No.	Colour	Code	Sales unit
33.	● 10 ch Black	190 1065 1	100
34.	○ 12 ch White	190 1265 1	100
35.	● 14 ch Green	190 1465 1	100
36.	● 16 ch Orange	190 1665 1	100

Catheter size – colour code system

CH/FR	Colour		Tube ^{mm}	Diameter
05CH/FR	Grey		1.7mm	○
06CH/FR	Light Green		2.0mm	○
07CH/FR	Violet		2.3mm	○
08CH/FR	Blue		2.7mm	○
10CH/FR	Black		3.3mm	○

Catheter size – colour code system

CH/FR	Colour		Tube ^{mm}	Diameter
12CH/FR	White		4.0mm	○
14CH/FR	Green		4.7mm	○
16CH/FR	Orange		5.3mm	○
18CH/FR	Red		6.0mm	○
20CH/FR	Yellow		6.7mm	○



COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company. © 2012 Covidien.

MS_ND_2012_EP13_0309

COVIDIEN (UK) COMMERCIAL LTD.
4500 PARKWAY,
WHITELEY,
FAREHAM, HAMPSHIRE,
PO15 7NY, UK

+ 44 (0) 203 027 1757 [T]
+ 44 (0) 203 027 1758 [F]
UKSALES@COVIDIEN.COM [E]

WWW.COVIDIEN.COM