

Dover™ Urine Drainage Bag 900010



Product Description	Dover™ Urine Drainage Bag
Bag capacity	2000ml
Product Packaging	
Items per case	1x20
Case dimension L, W & H	39.5 x 29.5 x 25 (cm)
No. of cases per pallet layer	8
No. of layers per pallet	8
Bar coding information case	GS1-128 barcode (01) 20884521177489
Product Materials	
Material bag	PVC
Material tubing	PVC 5401
Length of tubing	120cm
Material of inner packaging e.g. CSR wrap	70gsm paper, and 70 micron film
Sterilisation method	ETO
DEHP free	No
Latex free	Yes
Legal Requirements	
Country of manufacture	China
Shelf life	5 years
CE marked	Yes
Special Features	
NF sample port	Yes
Anti reflux valve	Yes
Air vent on drip chamber	Yes
Air vent on bag	Yes
Other	
Duration of use	24hrs (1 day) to 168hrs (7 days)