

MEDISANA®

Pulse oximeter PM 100

Art.-No. 79455 | EAN 40 15588 79455 1

Pulse oximeter to measure the saturation of oxygen in blood (SpO₂ %) and the heart rate (pulse)

View measurements on a OLED display

Display brightness adjustable over 10 levels

Simple „one touch“ operation – also ideal for travelling

OLED
DISPLAY



3 YEARS
WARRANTY

MDD 93/42/EWG



Application | Features

Pulse oximeter to measure the saturation of oxygen in blood (SpO₂ %) and the heart rate (pulse)

Oxygen saturation (SpO₂) and pulse rate can be displayed as bars or waves

View measurements on a OLED display

Viewable from 6 different angles

Display brightness adjustable over 10 levels

Simple „one touch“ operation – also ideal for travelling

Automatic switch-off after 8 seconds

Battery level indicator

Carrying strap included

Certified medical device*

Technical Data

Size: approx. 5,8 x 3,4 x 3,5 cm

Weight: approx. 53 g

Measuring range:
- SpO₂: 70% ~ 100%
- Pulse: 30 – 235 beats/min.

Batteries: 2 x 1,5 V (AAA)

Scope of Delivery

Pulse oximeter

Batteries

Instruction manual

Carrying strap

Logistical Data

Sales packaging

1 piece | tbd | tbd

Packaging unit

tbd | tbd | tbd

*Certification according to MDD (Medical Device Directive)