

AEROcase®



Emergency
Medical Service

Emergency bags | Emergency backpacks



Emergency Medical Service
Systems for professionals

AEROCASE®

Emergency bags | Emergency backpacks



AEROCASE® - Pro^{EMS} PXL1C Emergency backpack XL

Extremely durable emergency backpack. Designed for permanent daily use by emergency medical services. Provides safe transport of emergency equipment combined with maximum carrying convenience.

Spacious and practice-based interiors for a clearly arranged placement of substantial emergency equipment.

Well organized for professional emergency management.

2 Front pockets

- Interior fixing straps
- Plenty of space for large ampoule case, intubation kit, etc.

3 Side pockets

- For devices such as sphygmomanometer, stethoscope, sharps container, etc.

Document compartment

- For documents incl. clipboard

Processing details

- Zippers in signal coloring
- Wide hook-and-loop strips for patches
- Highly retroreflective reflective strips all around

The models

„Emergency backpacks XL, not just in size the greatest“

Made of AEROTEX® - Plan.

Measurements: 52 x 54 x 25 cm

Weight: 5.1 kg

Volume: approx. 55 L



Signal red



Luminous yellow



Equipment details

For professional use designed interior concept

Perfectly equipped for every kind of situation.

The modular and clearly arranged interior concept allows a safe and well organized placement of large emergency equipment:

- Hinged front and center walls with numerous useful fixing straps and compartments
- Structured, two-piece interior
- Integrated holder for:
0.8 L | 1 L | 2 L oxygen bottles
- Robust fixing strap (quality A):
 - Main compartment
 - Frame (encircling the inside)
 - Center compartment
 - Front pocket (bottom)



Comfortable transport

- Ergonomically formed handles with reinforced core
- Vertically and horizontally carryable
- Padded backstraps
- Storage compartment for backstraps
- Quick-fastening belts - fixation at stomach

Ideal Protection

- Hard-wearing AEROTEX® - Plan
- All around padded body
- Protected bottom area
- Trough shaped bottom
- 7 robust plastic pyramid feet



Individual organization

- 7 Module pockets in different sizes
- Transparent side and compartment
- With fastening eyes

AEROcase®



Emergency bags | Emergency backpacks

AEROcase® - Pro^{EMS} CL1 Emergency bag

Universal emergency bags for miscellaneous uses within emergency rescue services, doctors offices and other medical services.

Compact design accurately fitting:

- In trunk compartments of emergency vehicles
- In places with very limited space

The universal emergency bag concept

Processing details

- Handles and zippers in signal coloring
- Wide hook-and-loop strips for patches
- Highly retroreflective reflective strips all around

Related to practice and versatile

- Completely foldable (180 °)
- Ergonomical transport possibilities:
 - Formed handle with reinforced core
 - Optional carrying strap for comfortable shouldering
- Structured, practical interior layout with numerous options for organizing

The models

„Our emergency bag twins!“

For example used as: pediatric emergency bag / doctors emergency bag

Measurements: 45 x 34 x 19 cm

Weight: 2.9 kg

Volume: approx. 26 L



Signal red



Luminous yellow




Equipment details

Interior-Concept for the professional utilization


Customized organization made easy.

- Main compartment in the bottom with variable, removable separating strips
- Fixable, removable center wall with straps, loops and compartments in the upper main compartment for an ideal organization of equipment in fixed sizes
- Main compartment in the cover including lockable compartments and fixing straps (quality A)



Protected bottom and back areas

- Hard-wearing black plane material
- 4 solid plastic pyramid feet
- Protection from dirt, moisture and attrition
- Easy to clean and disinfect



Sealable

- Zipper comes with lockable sled

AEROcase® - Pro^{EMS} BOL1-B Oxygen bag

All around padded oxygen bag with storage room for:
O₂-bottles with pressure regulator | demand valve | resuscitation bag | O₂-masks and
further equipment for emergency resuscitation.

The answer for O₂-cylinders up to 5 liters.

Processing details

- Handles and zippers in signal coloring
- Wide hook-and-loop strips for patches
- Highly retroreflective reflective strips all around

The improvements

- Adapted coloring
- Concept for interior cleaning
- Easy-Access-Window
- Hook-and-loop strips for patches



*„Innovation meets tradition-
the classic with improved concept!“*

„Well proven and unmatched!“

Measurements: 55 x 21 x 31 cm

Weight: 1.4 kg

Volume: approx. 25 L



Luminous blue

Equipment details

Compact spatial wonder

- Large storage room for substantial oxygen equipment
- Main compartment with bottle holder and additional pockets
- Exterior front pocket
- Possibilities to hold: pressure regulator, demand valve, resuscitation bag, oxygen masks and further equipment for emergency resuscitation
- Clamping straps to restrict the max. opening of the bag
- Light gray plane material on inner side
- easy to clean and durable
- Additional front bag to house small items



Easy-Access-Window

- Permanent option to regulate the O₂-flow via easy-access-window with hook-and-loop system
- Continuous control of oxygen level

Convenient carrying

- Additional padded should strap
- Fastening straps for transport at a stretcher
- Hands free for more mobility



Inserted padding

- Protection mat of foam with 8 mm thickness



Integrated stable plastic plate

- For additional protection and stability



AEROcase®

Emergency bags | Emergency backpacks



AEROcase® - Pro^{EMS} BOXL1-B Oxygen bag XL

Ideal storage room | max. flexibility
max. comfort | Oxygen bag with holder for 2 L O₂-bottle
and all-in-one equipment concept.

The maximum for emergency resuscitation.

Processing details

- Handle and zippers in signal coloring
- Wide hook-and-loop strip for patches
- Highly retroreflective reflective strips

The improvements

- Adapted coloring
- Concept for interior cleaning
- Easy-Access-Window
- Hook-and-loop strips for patches



Example of filling

- **Main compartment:** O₂-bottle incl. pressure regulator | demand valve with 1.5 m hose and DIN connector | resuscitation bag | oxygen masks | intubation bag | etc.
- **Front bag:** Oxygen mask | device for NIV-CPAP ventilation | single use resuscitation device | etc.

„The all-rounder among oxygen bags!“

Measurements: 56 x 29 x 23 cm

Weight: 2.9 kg

Volume: approx. 26 L



Luminous blue

Equipment details

Flexible all-in-one solution

Sufficient space for substantial resuscitation equipment:

- Main compartment with holder for O₂-bottle via hook-and-loop strip
- Foldable front bag with hook-and-loop area and protective walls
- 4 holders on hook-and-loop basis for individual organization
- 3 parts with hook-and-loop areas:
Main compartment inside | Bottle holder | Front bag



Versatile carrying system

- Back straps in storage pocket on the backside
- Additional padded shoulder strap
- Fastening straps for transport at a stretcher
- Safe access when removing items

Easy-Access-Window

- Regulation of O₂-flow
- Control of oxygen level of the bottle



Ideal protection

- Hard-wearing AEROTex® - Plan
- All around padded body
- Protected bottom area
- Trough shaped bottom
- 3 robust plastic pyramid feet



Individual organization



O₂-bottle holder (plastic plate)



Loop holders



Storage compartment

AEROcase® - Pro^{EMS} CAS1 Small ampoule case |

CAL1 Large ampoule case

Ampoule cases for shatter-proof storage of ampoules, cannulas and syringes.

The ideal addition for bags and backpacks of the AEROcase® - Pro^{EMS} and AEROcase® - Pro^{1R} series.

Compact space saving solution

Small ampoule case

- Holds 24 ampoules (14 x 1-2 ml | 5 x 5 ml | 5 x 10 ml)
- Removable separating strip
- Integrated fixing loops on both sides
- Transparent pocket for list of contents
- All around integrated foam protects content from shattering
- Hook-and-loop strips on the backside for fastening of emergency pouches
- Space saving storage in emergency bags or backpacks



Spacious and safe

Large ampoule case

- Holds 60 ampoules
- Removable ampoule trays with professional ampoule strips and 3 elastic fixing loops
- Padded middle wall with spacious transparent compartments for syringes, cannulas, etc.
- All around integrated foam
- Hook-and-loop strips on the backside for fastening of emergency pouches
- Accurately fitting into the front bag of AEROcase® Emergency backpacks



The models

Measurements: 18 x 14 x 7.5 cm

Weight: 0.23 kg

Volume: approx. 1.8 L



Signal red

Measurements: 27 x 22 x 9.5 cm

Weight: 1.2 kg

Volume: approx. 5 L

AEROcase® - Pro^{EMS} CIS1 Intubation case

Intubation case made of high-quality and hard-wearing nylon that can be rolled up. The ideal equipment case for emergency bags and backpacks of the AEROcase® - Pro^{EMS} and AEROcase® - Pro^{1R} series.

Well organized emergency intubation.

High quality nylon

- Very tear-resistant fibers
- Hard-wearing and durable

Transparent pocket

- For a content list or something similar

Eye-catching quick-release buckle

- Secure lock
- Easy and quick opening

Eye-catching quick-release buckle

- Secure lock
- Easy and quick opening



2 Hook-and-loop strips

- Sichere Befestigung in Notfallbehältnissen

Well prepared intubation management

- Clearly arranged interior
- Solid, wide elastic straps (quality A)
- Spacious, separately closable inside pocket
- Secure and compact intake of: Laryngoscope, blades, tracheal tubes, etc.
- Flaps on the side protect the contents
- Foldable and therefore space saving



Belt adjustable in length

- Variable adjustments due to hook-and-loop fastener



Measurements: 30 x 16 x 4 cm

Weight: 0.24 kg

Volume: approx. 2 L

Luminous blue





Product overview Pro^{EMS} | Emergency bags and emergency backpacks


AEROCASE® - Pro^{EMS} PXL1C Emergency backpack XL

Item. No.	Description	Size	Material	Color	Dimensions (WxHxD)*	Weight	Volume
HT03-MPXL1C-R	AEROCASE® - Pro ^{EMS} PXL1C Emergency backpack 	XL	Plane	Red	52 x 54 x 25 cm	5.1 kg	approx. 55 L
HT03-MPXL1C-Y	AEROCASE® - Pro ^{EMS} PXL1C Emergency backpack 	XL	Plane	Yellow	52 x 54 x 25 cm	5.1 kg	approx. 55 L


AEROCASE® - Pro^{EMS} CL1 Emergency bag

Item. No.	Description	Size	Material	Color	Dimensions (WxHxD)*	Weight	Volume
HT03-MCL1-R	AEROCASE® - Pro ^{EMS} CL1 Emergency bag 	L	Plane	Red	45 x 34 x 19 cm	2.9 kg	approx. 26 L
HT03-MCL1-Y	AEROCASE® - Pro ^{EMS} CL1 Emergency bag 	L	Plane	Yellow	45 x 34 x 19 cm	2.9 kg	approx. 26 L


AEROCASE® - Pro^{EMS} B0L1-B Oxygen bag

Item. No.	Description	Size	Material	Color	Dimensions (WxHxD)*	Weight	Volume
HT03-MB0L1-B	AEROCASE® - Pro ^{EMS} B0L1, Oxygen bag 	L	Plane	Blue	55 x 21 x 31 cm	1.4 kg	approx. 25 L


AEROCASE® - Pro^{EMS} B0XL1-B Oxygen bag XL

Item. No.	Description	Size	Material	Color	Dimensions (WxHxD)*	Weight	Volumen
HT03-MB0XL1-B	AEROCASE® - Pro ^{EMS} B0XL1, Sauerstofftasche 	XL	Plane	Blue	56 x 29 x 23 cm	2.9 kg	approx. 27 L


AEROCASE® - Pro^{EMS} CAS1 Small ampoule case

Item. No.	Description	Size	Material	Color	Dimensions (WxHxD)*	Weight	Volumen
HT03-MCAS1	AEROCASE® - Pro ^{EMS} CAS1, Small ampoule case 	S	Plane	Red	18 x 14 x 7,5 cm	0.23 kg	approx. 1.8 L

AEROCASE® - Pro^{EMS} CAL1 Large ampoule case

Item. No.	Description	Size	Material	Color	Dimensions (WxHxD)*	Weight	Volumen
HT03-MCAL1	AEROCASE® - Pro ^{EMS} CAL1, Large ampoule case 	L	Plane	Red	27 x 22 x 9,5 cm	1.2 kg	approx. 5 L

AEROCASE® - Pro^{EMS} CIS1 Intubation case

Item. No.	Description	Size	Material	Color	Dimensions (WxHxD)*	Weight	Volumen
HT02-MCIS1-B	AEROCASE® - Pro ^{EMS} CIS1 Intubation case 	S	Nylon	Blue	30 x 16 x 4 cm	0.24 kg	approx. 2 L

*Dimensions including front and side pockets